

154

163

169

- Pince
- Pince
- Machine
- Pince
- Machine
- Pince

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SIGMA 1

COPYRIGHT KONTRON

All rights reserved. The information contained in this publication may not be used for any purpose other than that for which it was originally supplied. The publication may not be reproduced in part or in whole without written consent of Kontron. In order to maintain and improve standards of manufacturing, methods of functioning and to increase reliability, Kontron equipments are periodically reviewed. For this reason, the contents of this publication are subject to change without notice.

The equipment described is designed and manufactured by KONTRON S.A. (Montigny Le Bretonneux, France), a member of the Kontron group of companies.

QUALITY, RELIABILITY AND SAFETY.

This equipment has been designed with an emphasis on quality, reliability and safety, but Kontron can only accept responsibility for these aspects providing the following conditions are met :

- a) Electrical installations of the room or building in which the equipment is to be used must comply with regulations specified by the country in which the equipment is to be used.
- b) The equipment is used in accordance with the instructions for use provided by kontron(Operating manual).
- c) All modifications and repairs to the equipment are carried out by authorized Kontron personnel, their agents or authorized hospital technicians.

Your local Kontron company or agent is :

(To be entered by local Kontron company or agent.)

1. TABLE OF CONTENTS

1. TABLE OF CONTENTS	1
2. INTRODUCTION AND TECHNICAL SPECIFICATIONS	5
2.1. INTRODUCTION	5
2.2. CAUTION AND WARNINGS	5
2.3. CLEANING AND STERILISATION	6
2.4. TECHNICAL SPECIFICATIONS	7
3. THEORY OF OPERATION	13
3.1. BLOCK DIAGRAM	14
3.2. MAIN CLOCK AND SIGNALS CONTROL	17
3.4. SCON p.c.b.	18
3.5. MEMEX p.c.b.	24
3.6. PERIP p.c.b.	28
3.7. FIST p.c.b.	32
3.8. CARDIS p.c.b.	38
3.9. WOMOT p.c.b.	44
3.10. IFDOD p.c.b.	48
3.11. MEFRO p.c.b.	52
3.12. SISEND p.c.b.	56
3.13. TFOC p.c.b.	64
3.14. CABLING	68
3.15. VIDEO KIT	69
3.16. POWER SUPPLY	69
3.17. E.C.G. MODUL	69

4. ADJUSTMENT AND CALIBRATION (UNIT OPEN)	71
4.1. INTRODUCTION	71
4.2. CALIBRATION EQUIPMENT	72
4.3. TV MONITOR p.c.b.	82
5. SAFETY AND PERFORMANCE CHECKS (UNIT CLOSED)	91
5.1. TEST OF CLOSED MACHINE	91
5.2. SAFETY TEST	97
6. ASSEMBLY AND DISASSEMBLY	101
6.1. CMOS COMPONENTS	101
6.2. EXPLODED VIEW	101
7. FAULT FINDING AND STANDARD EXCHANGE	103
8. PARTS LIST	105
8.1. RECOMMENDED SPARE PARTS LIST	105
8.2. COMPLETE SPARE PARTS LIST	106
9. E.S.I.	
10. SCHEMATICS	
11. APPENDIX	

Introduction and Technical Specifications

2. INTRODUCTION TECHNICAL SPECIFICATIONS

2.1. Introduction

The service manual contains the necessary documentation for the servicing and repair of the SIGMA 1. It is assumed that readers are familiar with the operating manual inserted in chapter 0 of this document.

2.2. Cautions and warnings

NOTE: Maintenance may only be carried out by qualified service personnel.

The SIGMA 1 has CMOS parts. During maintenance it is essential to prevent transfer of static electrical charges to circuit boards and components by :

- grounding instrument
- grounding soldering iron
- wearing grounded wrist strap (via 1 MOhm to ground)
- adopting identical precaution for all CMOS spare components.

It is imperative to keep surface of the floating part of the E.C.G. p.b.c. between components and unit chassis insulated and absolutely clean.

The following general rules should be answered when operating or servicing SIGMA 1 :

- *The equipment must not be used in the presence of flammable anaesthetics.*
- *Operation of equipment in an environment with increased oxygen concentration is inadvisable.*

If the unit is opened it is imperative to perform all safety and functional checks described in chapter 5 of this manual after reassembly, to ensure that functional safety demands are met.

Functional checks after repair is required but not safety check, unless the power supply module, mains cabling, or the E.C.G. amplifier was repaired.

2.3. Cleaning and sterilization

2.3.1. CLEANING

Clean the surface of the Sigma 1 with a dry cloth. If more extensive cleaning is required, switch off the instrument and disconnect the power cord. Use a slightly dampened soft cloth and very mild detergent solution to clean exterior surfaces. Never use abrasive cleansers or steel wool, or extreme electrical disturbance may result.

With the Sigma 1 powered off, clean the transducer with a very mild detergent solution and a slightly dampened cloth.

If more extensive cleaning is required, contact your Service Support Representatives. Never remove the protective covers of the Kontron instrument. There are hazardous voltage levels inside.

2.3.2. STERILIZATION

Do not pressurize. Sterelize using non denatured alcohol and soft cloth. Never immerse the transducer in any liquid.

Sterilization may be carried out with paraformaldhehyde (polyoxymethylene) in tablets, at ambient temperature.

2.4. Technical specifications

2.4.1. COMMON SPECIFICATIONS FOR SECTOR AND LINEAR SCANNING

a) Image format

There are many different formats. They are shown in the enclosed table, which shows the individual combinations of imaging angles, depth, line density and frame rate. Note that it always takes an integer number of TV-frames to make one ultrasound image.

b) Ultrasound center frequency

Nominal frequencies : 3.5MHz, 5MHz.

Actual frequencies : very large band signal processing.

In the near range frequencies up to 9MHz are used with 3.5MHz sector transducer.

c) scale factor

On the built-in monitor, the image is always about 75 mm high. This results in the following scale factors :

imaging depth (cm)	scale factor (approximately)
6	1,25
8	0,9
10	0,75
13	0,6
18	0,4
23	0,3

Imaging depths between 8 and 23cm can be selected with 3.5MHz transducer. 5MHz and 7.5MHz transducers can be operated between 6 and 18 cm depth setting, although some among them may be useless because of limited acoustic penetration.

d) Gain

There are separate overall gain sliders for 2D image and for TM image (60 dB range).

The imaging depth is divided into 9 zones with 80 dB range in the first zone and 60dB range in the remaining zones.

e) Dynamic range

The dynamic range which is displayed on the screen can be controlled by the REJECT slider in the range from 25 dB to 60 dB.

f) Gamma-curve

The user can select from 4 different postprocessing curves. He can use normal video or inverse video.

g) Enhance

The user can select from 7 different types of enhancement (including "off").

Note : Gamma-curve and enhance functions are independently selectable for 2D and for TM

h) Filtering

The user can choose from 4 different filter functions trading signal to noise ratio and grey-scale against rendition of moving objects.

i) Scan converter features

The scan converter can display :

- 2D formats listed in a).
- Double 2D format. Each of the two images can be independently acquired using any of the formats available in the instrument. (Mixing of linear and sector images is also possible).
- Composite double 2D format in the case of linear scanning for displaying large objects.
- All formats combined with time motion.
(image is refreshed after one screen of TM).
- Full screen time motion mode at 10, 20 and 50 mm/sec.
- ECG trace on screen.
- ECG triggered image.
- Double ECG triggered image.
- Ciné mode for storing an image sequence.

j) Magnifier

The user can select a part of the image by the trackball and display it with a magnification factor of 2. This works with either real-time or frozen images.

k) Markers

There are cm-markers on the left side and on top of the image (respectively time markers in TM).

The position of the TM line is marked by a dotted line.

l) Alphanumerics on screen

- There is a full alphanumeric keyboard.
- There is a special field for the name of the patient.
- The user can write everywhere on the screen.
- He can move the cursor using the space/line feed keys, or using the trackball.
- The user can store a short text in the permanent memory and recall it with a keystroke.
- There is a technical data field which can be switched ON and OFF. It can display frequency, depth, near/far, high resolution/survey.
- There are four softkey fields at the bottom of the image, where the softkey functions are shown. The following features use the soft keys : Enhance, Filter, TM-speed, SET, ECG-sync, CALC.
- There is a symbol showing left-right inversion, freeze and magnification.

m) Grey scale

On the left side of the screen, there is a grey scale test pattern with 16 shades of grey.

n) Scan converter resolution

The visible portion of the screen has a resolution of 380 pixels horizontal and 512 pixels vertical with a resolution of 6 bits (60 shades of grey).

The graphics and alphanumerics display has a resolution of 380 x 256 pixels.

o) Set functions

There are 4 set functions :

- Program sound velocity used for measurements and display.
- Set clock circuit.
- Store short text into permanent memory.
- Store desired power up - state of machine (e.g. enhance and filter functions) into permanent memory.

p) Set state of machine

For linear and sector mode, there are 10 predefined machine states available in the permanent memory. They are activated by SET n where n is a number between 0 and 9. Four states are permanently set by the factory. Six states are user programmable.

q) Panic function

If the user is lost when manipulating signal processing modes, he can press SET I in order to return to the power-up state.

r) Measurements

Distance, circumference and area (in 2D).

Distance, time and slope (in TM).

Markers are moved by the trackball.

Double distance in 2D.

Two independent distances in double 2D.

s) Calculations

Built-in biometry functions are :

- for 0B/GYN : BPD, THD, ABD, FML, CRL, GES, HIP A, B.
- for cardiology : on TM → Distance ratio
Ventricular volume (Teichholz).
on 2D → Volume calculation
(SIMPSON or HEMI-ELLIPSOID).
- For internal medicine : RUV
In double 2D format, two independent calculations can be done, one on each image.
- User programmable : six functions based on distance, area or circumference input delivering a single output value of arbitrary dimension.

t) Clock

There is a permanent clock in the machine displaying the date.

u) Power consumption

Approximatively 200 W typical.

2.4.2. SECTORIAL PART PERFORMANCE

a) Image formats (50Hz line frequency)

SIGMA 1

FORMATS AND SCAN RATES - WOBBLER SCANHEADS (1)

Depth	3.5 MHz A						3.5 MHz C					
	HR			S			HR			S		
80	60	16 2/3	161	80	16 2/3	161	60	16 2/3	161	60	25	109
	90	12 1/2	203	80	16 2/3	167	80	12 1/2	209	80	16 2/3	167
	105	8 1/3	255	90	12 1/2	203	90	12 1/2	203	90	16 2/3	163
100	60	16 2/3	161	60	16 2/3	161	60	16 2/3	161	60	25	109
	90	12 1/2	203	80	16 2/3	167	80	12 1/2	209	80	16 2/3	167
	105	8 1/3	255	90	12 1/2	203	90	12 1/2	203	90	16 2/3	163
130	60	12 1/2	161	60	16 2/3	129	60	16 2/3	129	60	25	85
	90	12 1/2	163	80	16 2/3	127	80	12 1/2	167	80	16 2/3	127
	105	8 1/3	255	90	12 1/2	163	90	12 1/2	163	90	16 2/3	137
180	60	12 1/2	129	60	16 2/3	109	60	12 1/2	129	60	16 2/3	109
	90	10	163	80	12 1/2	127	80	12 1/2	127	80	16 2/3	107
	105	8 1/3	255	90	12 1/2	137	90	12 1/2	137	90	16 2/3	117
230	60	10	129	60	12 1/2	109	60	12 1/2	109	60	16 2/3	85
	90	8 1/3	163	80	12 1/2	107	80	10	127	80	12 1/2	107
	105	7	211	90	12 1/2	137	90	12 1/2	137	90	12 1/2	137

60 16 2/3 161
 angle frame rate number of lines

SIGMA 1

FORMATS AND SCAN RATES - WOBBLER SCANHEADS (2)

Depth	5.0 MHz B / 7.5 MHz B						5.0 MHz C					
	HR			S			HR			S		
60	60	16 2/3	161	60	16 2/3	161	60	16 2/3	161	60	25	109
	90	12 1/2	203	80	16 2/3	167	80	12 1/2	209	80	16 2/3	167
	105	8 1/3	255	90	12 1/2	203	90	12 1/2	203	90	16 2/3	163
80	60	16 2/3	161	60	16 2/3	161	60	16 2/3	161	60	25	109
	90	12 1/2	203	80	16 2/3	167	80	12 1/2	209	80	16 2/3	167
	105	8 1/3	255	90	12 1/2	203	90	12 1/2	203	90	16 2/3	163
100	60	16 2/3	161	60	16 2/3	161	60	16 2/3	161	60	25	109
	90	12 1/2	203	80	16 2/3	167	80	12 1/2	209	80	16 2/3	167
	105	8 1/3	255	90	12 1/2	203	90	12 1/2	203	90	16 2/3	163
130	60	12 1/2	161	60	16 2/3	129	60	16 2/3	129	60	25	85
	90	12 1/2	163	80	16 2/3	127	80	12 1/2	167	80	16 2/3	127
	105	8 1/3	255	90	12 1/2	163	90	12 1/2	163	90	16 2/3	137
180	60	12 1/2	129	60	16 2/3	109	60	12 1/2	129	60	16 2/3	109
	90	10	163	80	12 1/2	127	80	12 1/2	127	80	16 2/3	107
	105	8 1/3	255	90	12 1/2	137	90	12 1/2	137	90	16 2/3	117

60 16 2/3 161
 angle frame rate number of lines

b) Focussing

Static focus. A choice of different transducers is available with optimized focal distances for different applications.

2.4.3. LINEAR PART PERFORMANCE

a) Image formats

SIGMA 1

FORMATS AND SCAN RATES - ULAP TRANSDUCERS

Frame rate			Number of acoustical beams							
Depth	Scanhead									
	3.5 MHz Lin		5.0 MHz Lin		7.5 MHz Lin		3.5 MHz Cur1		5.0 MHz Cur1	
60 S			25	166	25	170			25	170
HR			16 2/3	249	16 2/3	255			16 2/3	255
80 S	16 2/3	196	25	166	25	170	25	170	25	170
HR	16 2/3	196	16 2/3	249	16 2/3	255	16 2/3	255	16 2/3	255
100 S	16 2/3	196	25	166	25	170	25	170	25	170
HR	16 2/3	196	16 2/3	249	16 2/3	255	16 2/3	255	16 2/3	255
130 S	16 2/3	196	16 2/3	166	16 2/3	170	16 2/3	170	16 2/3	170
HR	16 2/3	196	16 2/3	249	12 1/2	255	12 1/2	255	12 1/2	255
160 S	12 1/2	196	16 2/3	166	16 2/3	170	16 2/3	170	16 2/3	170
HR	12 1/2	196	10	249	10	255	10	255	10	255
230 S	12 1/2	196					12 1/2	170		
HR	12 1/2	196					8 1/3	255		
Nr. Seg.	108		93		96		96		96	
Scan-W.	107 mm		66 mm		51 mm		67 mm *		51 mm *	

* at skinlevel 48.4°

b) Focussing

Dynamic focus with 6 different focal zones at reception.
Near/Far focus selectable at transmission.

2.4.4. OPERATING PARAMETERS

- . Power source : 110V range : 90 ... 132 V_{eff}
220V range : 180 ... 264 V_{eff}
- . Temperature : 0 to 40°C
- . Humidity :
- . Warm-up time : 1 minute
- . Inputs and outputs : 1V_{pp} into 75
- . Video output to recorder : Z = 75
- . Video output to ext. monitor : 1V_{pp} into 75
- . Input for ECG electrodes (same as MICROMON)
- . 25 pole D-Connector for TM recorder and Doppler (specs of TM-output, see description of CARDIS board)
- . 25 pole D-Connector for optional serial interface
- . 37 pole D-Connector for optional parallel interface


Functional Description

3

3. THEORY OF OPERATION

3.1. Block diagrams

3.1.1. OVERALL BLOCK DIAGRAM

3.1.2. INTRODUCTION

SIGMA 1 is an ultrasound continued linear sector scanner for diagnostics in abdominal, obstetric and gynaecology applications.

Some of the boards from Sigma 1 SC had to be redesigned:

SCON 2 : A new Memory - Management. Unit was installed to receive more ROM-space for one microprocessor.

MEFROB : A new transmitter circuit was designed.

FIST : FIST was already prepared for nonlinear recursive filters. In Sigma I these filters are used, therefore FIST has to be manufactured with faster RAM's (120 ns selection) and two additional PROMs.

One board is completely new :

SISEND : This board contains the linear transmitter and receiver and is familiar to the ULAP – Board "SEND".
The differences are :

- New transmit - pattern logic,
- New hybrides,
- Low power dissipation.

One board is taken from ULAP:

TFOC

The other boards are taken from Sigma 1 SC:

PERIP
CARDIS
MEMEX
WOMOT
IFDOD

3.1.3. DESCRIPTION OF OVERALL BLOCK DIAGRAM

The SIGMA 1 was conceived aiming at low cost and modularity. SIGMA 1 was designed using the latest electronic circuit technology. It is highly modular and completely microprocessor controlled. Therefore, the system has a very good cost-performance ratio.

The heart of SIGMA 1 is the scan converter, which is built on just two boards. It converts the incoming ultrasound echo information from any scanning format into a standard television format. It includes two overlay memories and such features as magnifier and time motion on screen even combined with a 2D-image. It includes two 8 bit microprocessor systems, one for real time control functions and the other for keyboard and overlay memory management. The compact design of the Scan converter was only possible by the application of "Gate Array" technology. Two semicustom integrated circuits were developed, each representing the functional equivalent of about 80 standard IC-packages. This offers important advantages in space, power consumption, reliability and production cost.

The basic building block of the analog sectorial part of SIGMA 1 is the IFDOD board. On it, there are the largest signal processing circuits like dynamic filters, logarithmic amplifier, dynamic range control and various "enhance" functions.

To that, we must add the motor control board (WOMOT) and the Transmit/Receive board (MEFRO) in order to get a mechanical sector scanner.

The basic building block of the analog linear part is the SISEND board which generates the optimized signals to drive the linear transducer elements for transmission and connects the signal of the required phase onto I/O lines to the ODD and EVEN boards in the transducer. It provides the initial amplification and phase selection of the received signals for transmission to the TFOC board.


The FOC board processes the received signals for apodization and focus.

The image quality of the scanner has been further improved by adding the frame filter board (FIST).

The gain control voltage is generated on the CARDIS board. On the same board, there are the circuits required for cardiology features : the interface to the ECG amplifier module and the TM recorder output..

3.2. Main clock and signals control

3.2.1. MAIN CLOCK AND SIGNALS CONTROL BLOCK DIAGRAM


3.2.2. DESCRIPTION OF MAIN CLOCK AND SIGNALS CONTROL BLOCK DIAGRAM

The whole system is driven by a 16 MHz master clock on SCON. For the bus the clock is divided by two. The resulting 8 MHz clock is called MCLK and is used by most boards as a time base.

The scan of an image is defined by two signals, which in case of a sector scan are produced by WOMOT: The ZBM signal is active during the whole image. The SRV signal triggers the individual transmit-pulses, based on the actual position of the mechanical scanhead. In case of TM or electronic linear scanning ZBM is not needed anymore. WOMOT is disconnected from the bus and SRV is driven by a pulse called BCLK (coming from RGA on SCON). The gain control voltage is delivered by the CARDIS board through the GCVext line.

3.4. Scon p.b.c.

3.4.1. SCON BLOCK DIAGRAM

3.4.2. DESCRIPTION OF SCON 2

a) SCON-Scan converter

The SCON2 board together with the MEMEX board is a complete controller/scan converter. There are two microprocessors : one for real time control of the scanning and the other for keyboard and overlay memory management. Video synchronization and the read and write addresses of the image memory are generated by two "Gate Arrays" (semicustom integrated circuits). Each of them represents an entire board of TTL-logic!

At the signal input, there is a 6 bit flash A/D converter which is not inserted : it is only required, if there is no FIST board in the system. Else, the digital output of FIST is used directly. Then, the signal goes through the "Data PROM", which serves as multiplexer for image initialization and as look-up table in order to do some processing in TM-mode. Subsequently, it is fed into the image memory on MEMEX. The output from the image memory goes to the interpolation stage. It detects not valid pixels. Then, the signal goes through the "Postprocessing PROM", where different gamma-curves are stored (including video inversion). The PROM output has 8 bit-resolution and is converted by a 8 bit digital to analog converter. The resulting analog signal is smoothed by an analog filter, which improves the interpolation performance and the appearance of magnified images. Then, the signal is mixed with video SYNC signals and overlay, and blanking information is added. Finally, there are two 75 Ohm output stages.

b) Microprocessor N°1


This is a 8085 based microprocessor dedicated to all real time control tasks and some image initialization tasks. For instance, before scanning it must fill the zone of a sector image by the "not valid" code in order to enable the interpolator to detect pixels, which were not written by the scanning.

The processor has a bus connection to the Analog Part and the Frame Filter. The address-space organization of this connection is given in *Table 1* (following page). The organization of the whole address space is shown in table SCON.MEMORY SPACE 2, but because address decoders are PALs, this may change during software development. The processor controls the "Read Gate Array" and the "Write Gate Array".

The restart inputs are used for synchronizing to the television retrace RTR (RST 6.5), to the Send/Receive clock SRV (RST 7.5) and to the frame start/stop of the mechanical scanhead ZBM (RST 5.5).

The microprocessor N°1 gets its basic information from the microprocessor N°2, which can put the microprocessor N°1 into the HOLD state and overtake its bus. The microprocessor N°2 does so usually during the vertical television retrace, on condition that microprocessor N°1 is not scanning an image. This overtake requires some handshakes between microprocessor N°1 and microprocessor N°2. The microprocessor N°1 gets the RTR interrupt, prepares the microprocessor N°2 bus access and interrupts microprocessor N°2 by setting the SOD output to a high level. Microprocessor N°2 gets that signal on RST 5.5 and overtakes the bus of microprocessor N°1 until the end of the retrace is given by a RST 7.5 on microprocessor N°2. Part of the RAM space of microprocessor N°1 is used as a mailbox (PATCH-PROM), in order to exchange information between microprocessor N°2 and microprocessor N°1.

TABLE 1


FIGURE : SCON MEMORY-SPACE.2 I/O $\mu P N^{\circ}1$

c) Microprocessor N°2

This is a 8085 based microprocessor, which performs all back-ground tasks of the system. On SCON2, there is a triple timer (used e.g. for the "beep"). Microprocessor N°2 has a bus connection to MEMEX for writing into the overlay II. It has a reduced bus connection (256 bytes) to PERIP, where the 8-channel A/D converter, two keyboard controllers, RS 232 and Floppy interfaces and a real time clock/permanent RAM are located. *Table 2* SCON.MEMORY Space 3 gives the memory space organization. It is also decoded by programmable devices and may therefore be changed during software development. The organization of the I/O space can be used from the description of PERIP (see table PERIP.CONTROL).

The main software functions of microprocessor N°2 are keyboard and gain control management, ECG tracing and trigger and overlay memory control. Because ROM-space of microprocessor N°2 is too small for all required functions, it can use 16k from the microprocessor N°1 space (PATCH-PROM) through the bus coupler for programs, Which are only used during frozen image. In this case, microprocessor N°1 is set to the HOLD state.

TABLE 2

- . The 16K byte PATCH PROM can be replaced by a 8K byte static RAM for RAM intensive applications.
- . An additional ROM-space of maximal 56K Byte is placed on the memory space of Microprocessor N°2 and therefore accessible at any time.

It is divided into 4K Byte-Blocks (programmable in a PAL) which are selected by the same window-selector-bytes in "Status 2" as for the PATCH-PROM, setting Bit 3 to 1 instead of 0.

The devices are selected by jumpers ST 1 to ST 4. The setting of these jumpers is shown in *Table SCON JUMPERS*.

There are 8 kbytes of static RAM with built-in battery backup. This nonvolatile memory is used for storing the user biometry functions and the state setups.

The memory space of SCON2 is shown in *table 2* page 21. The meaning of the bits of the SCON-Status-Registers (one on Microprocessor N°1 - side, one on N°2 - side) is shown in *table SCON Status Registers* (following page).

Finally SCON2 is prepared for use with a 8MHz processor-clock. In this case ST 5 (see *table SCON JUMPERS*) has to be programmed to include WAIT-cycles for all peripheral accesses and PERIP has to be fast enough to match the speed.

However the microprocessor is still running at 5.33MHz in Sigma 1.

ST 1	: Select the device which is installed at the middle socket on microprocessor 1 side (Z 110) :
	1 - 2 : RAM 8 K x 8 static
	2 - 3 : EPROM 16 K x 8
ST 2 to ST 4 :	Select the size of EPROM installed in the sockets Z 512, Z 513, Z514 on microprocessor 2 side :
	1 - 2 : EPROM 16 K x 8
	2 - 3 : EPROM 32 K x 8
ST 5	: Determines the microprocessor-clock for both microprocessors
	"8" : 8MHz clock (results in 250 ns cycle-time)
	"5" : 5.33MHz clock (results in 375 ns cycle time)

SCON JUMPERS Setting of Jumpers table

SCON STATUS micro-Processor 1 : Bit 0

- 1 DATA PROM
- 2
- 3 WE (OV)
- 4 WE (R)
- 5 INPOLEN
- 6 Ping Pong Status Plane A 0 = Write
- 7 Ping Pong Status Plane B 1 = Read

SCON STATUS micro-Processor 2 : Bit 0


POST PROCESSING

- 1
- 2 n.a.
- 3 MMU-Control : 0 = micro Processor
N°1 Patch
1 = micro Processor
N°2 Patch
- 4
- 5 Window pointer (for micro processor
N°1-side)
- 6 Patch selector (for micro processor
N°2-side)

SCON STATUS REGISTER

3.5. Memex p.c.b.

3.5.1. MEMEX BLOCK DIAGRAM

3.5.2. DESCRIPTION OF MEMEX BLOCK DIAGRAM

Memex is an extension of the SCON board. It is directly connected to it by two 64 pin connectors. There are three different screen memories on SIGMA 1.

The largest one is the image meory which consist of two blocks of $512 \times 512 \times 6$ bits each. This is due to the "PING-PONG" scan conversion principle, where one memory block is read according to TV-standard, while a new image is written into the other block including all required geometric corrections. Then, the blocks are exchanged. The SCON board supplies two sets of address and control lines for read and write respectively. Each memory block has a multiplexer which selects either read or write address depending on the Ping-Pong status bits.

The second memory is the OVERLAY I, which has $512 \times 256 \times 2$ bits capacity. With two data bits, we can encode transparent, black, white and flashing pixels.

This memory is permanently connected to the read address lines. Therefore, the readout occurs in parallel with the image memory. The microprocessor can only access this memory during the vertical TV-retrace. The access-circuit is part of the "Read Gate Array".

The third memory is the OVERLAY II, which has $512 \times 256 \times 1$ bit capacity. It is far better accessible than OVERLAY I, because the microprocessor can access it during the whole TV-frame (but not during TV-retrace). After all, the access is not fully transparent, but produces transparent pixels on the TV-ouput. Therefore, flickering can occur, if many white pixels are displayed on OV II.

This memory is mainly used for displaying the ECG-trace and the TM-cursor, because this requires fast memory acces. Alphanumerics, centimeter scales and calipers are displayed by OVERLAY I. The microprocessor can not read from OVERLAY II. The opportunity to write into the image memory by the microprocessor is used for generating the greyscale testpattern.

Figure 3 shows how the image memory is organized in order to have a readout rate of 8MHz and sampling clock up to 5.33 MHz, using dynamic RAMs with 150 ns access and 270 ns cycle time. During read, a 2-word PAGE MODE is used and two memory banks are interleaved. During write, two words are written in parallel into two memory banks. Figure 4 shows the timing diagrams for read and write.


FIGURE 4




IMAGE MEMORY READ TIMING

Note: The frequency of SCK is variable (depending on display size, scale and angle). The first SCK cycle of 187ns is shown.

3.6. Perip p.c.b.

3.6.1. PERIP BLOCK DIAGRAM

3.6.2. DESCRIPTION OF PERIP BLOCK DIAGRAM

The watch dog detects the rare occasions where the microprocessor does not run properly. Then, it makes an automatic reset. There is a battery powered CMOS circuit containing a clock and a permanent memory of 56 bytes.

The board distributes the voltages from the power supply to different locations. The digital part is powered through a cable to the SCON board. The analog part is powered by the digital part. PERIP provides also the supply connection for the fans by a temperature dependent inside of SIGMA ISC. This minimizes fan noise. The microprocessor bus is interfaced to the keyboard and footswitches by two keyboards controllers. The trackball is connected via a trackball-interface, which is built by two PALs (field programmable array logic). In this interface, the clock pulses which are generated by the trackball as a function of the ball-rotation are counted. After several counts, an interrupt request is sent to the microprocessor, which can read the state of the counter. Some control signals which are required for the TM/Doppler interface, are driven or read by a I/O Port. These signals are listed in *table 5*.

There is an eight channel A/D converter. It makes it possible to access by the microprocessor the ECG-voltage and several slider or potentiometer voltages. Finally, there is a room for fixing a small piggy-back circuit on PERIP called PARSIO (for "parallel/serial input/output"). It contains a parallel and a RS 232 interface. The control words of PERIP are listed in *Table 6*.

TABLE 5

ADDRESS (HEX)	FUNCTION
00 - 3F	REAL TIME CLOCK WITH MEMORY
40 - 4F	KEYBOARD CONTROLLER I
49 - 4F	KEYBOARD CONTROLLER II
50 - 58	ANALOG TO DIGITAL CONVERTER
59 - 5F	STATUS BYTE
60 - 68	RS 232 INTERFACE
69 - 6F	PARALLEL INTERFACE
70 - 78	TRACKBALL INTERFACE
79 - 7F	NOT USED


TABLE 6

ADDRESS (HEX)	BIT	FUNCTION
59 (STATUS WRITE)	0	PAGE PRINT COMMAND (ACTIVE HIGH)
	1	BLANK
	2	TRACKBALL INTERRUPT ENABLE
	3	GATE ENABLE (ACTIVE LOW EXT MARKS)
	4	FBEN
	5	ECG TRIG INTERRUPT ENABLE
	6	PARALLEL INTERFACE INTERRUPT ENABLE
59 (STATUS READ)	7	PARALLEL INTERFACE RESET
	0	ECG INTERRUPT PENDING
	1	ECG FAULT
	2	STRAP STATUS PAL/NTSC
	3	CARDIS (DETECTS PRESENCE OF BOARD)
	4	SCAD
	5	TRAD
59 (STATUS READ)	6	O
	7	PARALLEL INTERFACE INTERRUPT PENDING

ADDRESS (HEX)	A/D CONVERTER SIGNALS
50	FILTERED ECG
51	GAIN IMAGE
52	GAIN TM
53	SLIDER I
54	SLIDER X
55	NOT USED
56	0.67*FAN SUPPLY VOLTAGE
57	REJECT

3.7. Fist p.c.b.

3.7.1. FIST BLOCK DIAGRAM

3.7.2. DESCRIPTION OF FIST BLOCK DIAGRAM

This board contains two digital filter circuits. Its purpose is to increase image quality and to serve as a buffer memory for special functions like double image and Time Motion.

The first filter is a "beam filter". It can be switched ON and OFF by microprocessor-control. The beam filter calculates the linear average of two consecutive ultrasound beams. This reduces image granularity and increases signal noise ratio. The circuit contains $1k \times 8$ fast static RAM for the storage of the information of one ultrasound beam.

The second filter is a "frame filter". It can perform various filter functions (selected by the microprocessor) by mixing the echo amplitude of consecutive ultrasound images. Therefore, it contains a frame memory for images with up to 256 beams of 512 samples each.

The frame filters, which are actually implemented on the board, are defined by the PROM contents. There is a choice of "averaging" filters which gradually improve the smoothness of the image but also tend to smear moving objects.

In addition, there is one MTI filter function. At the input of the board there is a 6 bits analog to digital converter. It converts the incoming echo signal (called LFOUT, coming from the IFDOD board) into a digital number. The reference voltage of the A/D converter must be trimmed to 3.0 Volts. The sample clock (SCK) is delivered by the SCON board. Its frequency is variable with imaging scale and angle. The minimum sample interval is 187 ns.

For the synchronization of the filter, it is important to know the interval, during which the beam echoes are written into the image memory on SCON. For this purpose the image-memory write-clock CAS W is transmitted from SCON to FIST.

The pixel addressing is generated by a 12 bits counter on FIST by counting SCK-pulses. The beam addresses for the frame memory are controlled by the real-time microprocessor of SCON. The short data cycle of 187 ns required the interleaving of two banks of dynamics memories for the frame filter. For the definition of the controlwords, see *Tables 7 and 8* (following pages). Timing diagrams are given in *FIST 2.* and *FIST 3.*


TABLE 7

ADDRESS	BIT	FUNCTION			
00	0	-			
	1	-			
	2	-			
	3	COL 0			
	4	.			
	5	.			
	6	.			
	7	.			
01	0	.			
	1	.			
	2	.			
	3	COL 9			
	4	-			
	5	-			
	6	-			
	7	-			
02	0	P0	0: BEAM FILTER OFF	1: BEAM F. ON	
	1	P1	00: FREEZE	01: FRAME F. OFF	
	2	P2	10: FRAME F.ON	11: NOT USED	
	3	P3	0: NORMAL	1: STORE FOR FREEZE	
	4	P4			
			0...3: PROM PAGE SELECT		
	5	P5			
	6	P6	00: RECURSIVE PROM	01: RECURSIVE x 0.5	
	7	P7	10: ZERO (TRANSVERSAL)	11: RECURSIVE x 0.5	
03	0				
	2				
	3				
	4	CONTROL WORD FOR MICROPROCESSOR INITIALIZED WITH 80H			
	5				
	6				
	7				

TABLE 8

VALUE (HEX)	FUNCTION OF FRAME FILTER	
0A	FILTER OFF (INCLUDING STORE FOR FREEZE)	
04		0
14	PROM-RECURSIVE	1
24	PROM PAGE	2
35		3
44		0
54	LINEAR RECURSIVE	1
64	PROM PAGE	2
74		3
84		0
94	TRANSVERSAL	1
A4	PROM PAGE	2
B4		3
08	LOGICAL "OR" OF ONE OF THE VALUES GIVEN ABOVE WITH 08 GIVES "STORE FOR FREEZE"	
01	THE VALUES GIVEN ABOVE RESULT IN "BEAM FILTER OFF". LOGICAL "OR" WITH 01 GIVES "BEAM FILTER ON".	
00	FRAME FILTER FREEZE	
NOTE : BEFORE ENTERING "FREEZE", ONE (AND ONLY ONE) FRAME MUST BE MADE IN THE STORE FOR "FREEZE" MODE.		


FIST 2

Note: Because of interleaved control signals, Dout A and Dout B are never active at the same time. Therefore, they can be directly connected.

Memories work in a read-modify-write cycle of 375ns min.

FIST 3


CRITICAL TIMING :

A: data hold at latch input after $\overline{\text{HOLD}}$	(min 20ns	LS 373)
B: write data setup before $\overline{\text{WE}}$	(min 20ns	6148-55)
C: address hold after $\overline{\text{WE}}$	(min 5ns	6148-55)

NOTE: For beam-filter off, $\overline{\text{CS}}$ of RAM is high and $\overline{\text{ENI}}$ is always low.

3.8. Cardis p.c.b.

3.8.1. CARDIS BLOCK DIAGRAM

3.8.2. DESCRIPTION OF CARDIS BLOCK DIAGRAM


The CARDIS board contains all functions which are required for a scanner which is specialized for cardiology. Basically, there are three sections : gain control, TM-recorder interface and ECG interface. The gain control circuit works with nine rectilinear potentiometers for time gain compensation. The zones are automatically distributed over the selected imaging depth by the same programmable control voltage that controls the TM deflection signal. There are two sliders for independent settings of overall gain in 2D and TM-mode. For optimum user comfort, there are two special features in the TGC-control. First the gain control range of the first slider is larger (80 dB) than the range of remaining sliders which is 60 dB. Second the settings of sliders 3 to 9 are not connected directly to the gain control circuit. Instead the voltage of two adjacent sliders are averaged by two resistors. This gives a smoothing effect on the TGC-curve which makes it easier to set a good curve. These two features are realized by resistor networks sitting on the POLIN board. POLIN is a small PCB behind the front panel, which carries the sliders.

The TGC-voltage is generated by an integrator. At the end of a reception, the integrator is preset to a value given by one of the overall gain potentiometers.

During the SRV-pulse (which is the preparation time for the next ultrasound beam) "DEPTH MUX" selects the voltage of slider 1 and integrates it over one depth interval in order to get the start gain. During reception, the integrator integrates during 8 intervals the voltage difference between two adjacent slider signals. The "DEPTH MUX" steps through sliders 1 to 9. It is controlled by a counter which in turn is clocked by a VCO. The frequency of this VCO is set accordingly to the selected penetration depth.

Figure 9 shows the timing diagrams of the TGC-circuit. The TM-interface provides two signals : TMZ and TMX. The TMZ signal is the echosignal which comes on LFOUT (driven by IFDOD) amplified and mixed with blanking and time-depth markers and the ECG channel signal LF EXT. This will be used for displaying the results of an external Doppler system on the screen of SIGMA 1.

FIGURE 9

The signal amplifier includes a high frequency enhancement which gives better TM quality. The output amplifier has strap selectable levels (see *Table 10*). The TMX-signal is the deflection ramp. It is generated by an integrator. Its slope is programmable depending on the selected penetration depth. The retrace has a constant slope of several hundred microseconds. The output amplifier also has strap selectable levels.

TABLE 10

STRAP 1 : TM-Z (SIGNAL)			STRAP 2 : TM-X (DEFLECTION)	
BLANK	WHITE	BLACK	DEFLECTION VOLTAGE (V_{pp})	
A-0.3	0	3	A	3
B 0.4	1	2.5	B	6
			C	16

Note : SIGNAL-Voltage measured at 50 Ohm termination. For best TM-quality, it is recommended to terminate the signal with 50 Ohm either inside of the TM-recorder or by a termination resistor at its input.

The time-depth markers are controlled by a programmable timer. The markers are tiny black dots surrounded by white rectangles. Horizontal dot size is given by a PAL where two modes are selectable by the microprocessor. The length of the dot pulses is automatically adjusted depending on the deflection speed in order to give approximately equal dot width for all depths.

Figure 11 shows the TMZ and TMX-signal including time-depth marker.

FIGURE 11



Figure 12 shows the example of a TM-output, showing tic-marks every 200 ms and time-depth marks every second. Note that the TM-interface must work at a 1KHz pulse repetition frequency on SRV. The ECG signal is preamplified and isolated by the ECG module, which is manufactured by Kontron Instruments Ltd, England. It is a modified version of the module used in MICROMON.

FIGURE 12

The module needs a 72 kHz clock and $\pm 10V$ power, which is supplied by CARDIS. The output signal of the module needs some filtering. On CARDIS, there are two notch filters at 50 Hz and at 60 Hz in order to remove interference from the mains voltage. There is a low pass with 20 Hz cutoff-frequency. Afterwards, the signal is amplified with a programmable gain by a multiplying D/A converter. The output signal is sent at four places. First to the ECG trigger section, second to the circuit which displays it on the TM-record, third to the A/D converter located on PERIP in order to display it on the TV-screen and finally to TM/Doppler connectors on the rear panel in order to display it by external means. Most of the signal on that connector come from the CARDIS board. The pinout of it is given on Table 14. The microprocessor control words are explained in Table CARDIS.CONTROL 13.

TABLE 13

ADDRESS (HEX)	BIT	FUNCTION
20	0	TIMER 0 (DIVIDES THE 1kHz SRV BY TEN, 100 Hz CLOCK)
	1	
	2	
	3	
	4	
	5	
	6	
	7	
21	0	TIMER 1 (NOT USED)
	1	
	2	
	3	
	4	
	5	
	6	
	7	
22	0	TIMER 2 (TRIGGERS DEPTH MARKERS)
	1	
	2	
	3	
	4	
	5	
	6	
	7	
23	0	CONTROL WORD OF TIMER
	1	
	2	
	3	
	4	
	5	
	6	
	7	
24	0	SLOPE OF TM-X DEFLECTION
	1	
	2	
	3	
	4	
	5	
	6	
	7	

TABLE 13


30	0	GCV DRIVE (ACTIV HIGH)
	1	TM-MODE (ACTIVE HIGH)
	2	DRIVE LFOUT BY EXTERNAL LF IN (ACTIV LOW)
	3	FORCE SLIDER SUPPLY TO 10 VOLTS REF (ACTIV LOW)
	4	DRIVE SRV BY EXTERNAL INPUT (ACTIV LOW)
	5	ENABLE ON-BOARD SRV (ACTIV HIGH)
	6	MARKER MODE SELECT
	7	TM-PROCESSING (ACTIVE HIGH)
34	0	
	1	
	2	
	3	
	4	ECG GAIN
	5	
	6	
	7	

TABLE 14

PIN	SIGNAL	FUNCTION
1	MCLK	8 MHz Master Clock output
2	14 GND	
2	BLCK	Beam Clock output
3	15 SRV int	Internal SRV output
3	GND	
4	16 SRV ext	External SRV input
4	ZMB	Zero Beam output
5	17 GND	
5	SCAD	Scan disable inputs to
6	18 TRAD	Transmit disable microprocessor
6	+5	Protected 5 Volts supply 0.3A max
7	19 FBEN	Fixed beam enable output
7	GATE	Gate-marker input
8	20 GND	
8	LF IN	External LF input
9	21 GND	
9	TMX	Line scan recorder X-signal output
10	22 GND	
10	TMZ	Line scan recorder Z-signal output
11	23 GND	
11	GCV	Gain control voltage input/output
12	24 EGC	ECG signal output
12	PAGE PRINT	Line recorder PRINT command
13	25 GND	
13	BLANK	Software controlled blank

3.9. Womot p.c.b.

3.9.1. WOMOT BLOCK DIAGRAM

3.9.2. DESCRIPTION OF WOMOT

The WOMOT board controls the servo system of the mechanically wobbling sector scanhead. There is a programmable function generator which generates the desired position control voltage. There is a servo loop, which keeps the error between the desired position and the actual position (which is measured by a position-sensor in the scanhead) to a minimum.

The trigger pulses for the acoustic transmission are derived from the actual position. For that purpose the latter is converted to a 10 bit digital value which addresses a PROM look-up table, where the transmit-positions are stored. There are security circuits which protect against double triggering and excessively short intervals between transmit pulses.


There is also a microprocessor input for reading the transducer code.

The function generator is driven by a programmable clock frequency in order to select the desired frame rate. A Modulo 3200 counter counts this clock and addresses a PROM table, where a sine-function is stored with 12 bit precision. The output of the PROM is converted into an analog voltage. The amplitude of the sinewave is adjusted by applying a programmable voltage to the reference input of the D/A converter. The phase of the sinewave can be controlled by the microprocessor by presetting the counter to a programmable state. To this sinewave, an offset voltage is added, which is also programmable by the microprocessor.

In the servo loop, the position sensor output is first multiplied by a scale factor voltage, which is delivered by the scanhead. The resulting voltage is fed to the A/D converter. The signal is also low pass filtered and subtracted from the function generator output. The resulting error is amplified by a proportional and a differential channel and subsequently drives the power servo amplifier. The latter works as current source, because this avoids any phaseshift due to the motor inductance. A security circuit controls the power amplifier. If an illegal state is detected, the motor current is interrupted. This protects the power amplifier and the scanhead from overheating.

Figure 15 shows the principle of the wobbler scanhead.

FIGURE 15

Tables 16 and 17 explain the microprocessor control words.

TABLE 16


ADDRESS (HEX)	BIT	FUNCTION
04	0	LOAD COUNTER PHASE
	1	
	2	
	3	
	4	
	5	
	6	
	7	
05	0	INPUT TRANSDUCER CODE
	2	
	3	
	4	
	5	STATE OF POWER UP RESET (1 = ACTIV)
	6	
	7	
07	0	DRIVE SRV LINE (ACTIV HIGH)
	1	RESET SECURITY CIRCUIT (ACTIV HIGH)
	2	SINE-LOOK-UP PROM: A12
	3	POSITION VOLTAGE LOW PASS CONTROL : NORMALLY LOW
	4	SCAN CASE (STATE PROM)
	5	
	6	
	7	
22	0	LOW BYTE
	1	
	2	
	3	
	4	
	5	
	6	
	7	
23	0	POSITION OFFSET CONTROL
	1	HIGH BYTE
	2	
	3	
	4	
	5	
	6	
	7	

TABLE 17

ADDRESS (HEX)	BIT	FUNCTION
24	0	TIMER 0 : SRV PULSE LENGTH
	1	
	2	
	3	
	4	
	5	
	6	
	7	
25	0	TIMER 1 : RECEIVE-TIME PROTECTION
	1	
	2	
	3	
	4	
	5	
	6	
	7	
26	0	TIMER 2 : PRESCALER FOR FUNCTION GENERATOR CLOCK
	1	
	2	
	3	
	4	
	5	
	6	
	7	
27	0	TIMER MODE CONTROL
	1	
	2	
	3	
	4	
	5	
	6	
	7	

3.10. Ifdod p.c.b.

3.10.1. IFDOD BLOCK DIAGRAM


3.10.2. DESCRIPTION OF IFDOD BLOCK DIAGRAM


IFDOD is the heart of the analog part. There is an interface between the microprocessor control bus and the frontend boards. It desactivates all digital control lines inside the analog part all the time, except when an access to the analog part occurs. This access must not occur during the reception time. That way, the risk of coupling digital noise into the receiver is minimized. Second, of all, there is a programmable gain control voltage generator for non-cardiac applications. It has three control parameters : near gain, far gain and overall gain. The change over between near and far gain is given by the hardware and has two selectable states for a coarse adaptation to the penetration depth. On the output, there is a multiplexer which can alternatively select an external gain control voltage. CARDIS delivers on this line its special gain control voltage. Based on the clinical evaluation of Sigma 1 it was decided, not to use the gain control circuit IFDOD. Hence for all application, the TGC curve is delivered by CARDIS. The circuit on IFDOD is only used for generating the TM-greyscale. In this cas, the "gain control voltage" is injected to the LFout line by a multiplexer which is located in front of the enhance circuit. Third of all, there is some logic for the real time control signals. There is a driver, which can drive the SRV-line by the Beam clock signal. There is also an EXOR gate, which generates microprocessor interrupts on the ZBM-transmissions.

The most important circuit on the board is the signal processing chain. First, the high frequency echo signal coming from the swept gain amplifier is filtered, in order to optimize the signal to noise ratio. There is a dynamic bandpass filter. It is explained in *figure 18*. In the near range, where the echo signals are strong and contain a lot of high frequency components, the center frequency of the filter is high and its bandwidth is large. With increasing depth, the high frequency components are strongly attenuated by the tissue. Therefore, it is optimal, to shift the center frequency of the filter to lower values and decrease the bandwidth. The voltage ramp which controls the filter is programmable and there are three different filters. This way, the filter function can be perfectly adapted to the transducer frequency and the desired penetration depth. After filtering the signal goes to the logarithmic amplifier. It is optimized for 60 dB dynamic range and clean impulse response from 2 to 10 MHz operating frequency.

FIGURE 18

QUALITATIVE CHARACTERISTIC OF DYNAMIC FILTER

After detection, the signal is low pass filtered. There are two selectable cutoff-frequencies. The next stage offers a choice of eight different gain and offset settings. This is used for the "REJECT" function and defines the dynamic range which is used for further processing and display. Now, the signal is processed by the "sharpening" or "enhance" circuit. Its purpose is to enhance the visibility of transitions in the image and eventually perform an automatic gain control. The optimal sharpening depends on the application of the instrument and on the user. Therefore, the user can select from seven states. There is an "off" state. There are three "expand" states, which enhance edges on the image without touching the grey levels of larger regions. This is suitable for cardiac imaging. There are three "compress" states, which produce sharp transitions but also decrease the contrast of large regions and include an AGC function. This is suitable for abdominal imaging.

The "expand" states are indicated on the screen as ENH1 to ENH3. The "compress" states as ENH-1 to ENH-3.

The output amplifier drives the LFOUT line. There is a multiplexer which can apply the gain control voltage to LFOUT. This is useful for generating a grey scale on the LFOUT line e.g. for adjusting the TM recorder.

The microprocessor control-words are given in *Table 19*.

TABLE 19

ADDRESS (HEX)	BIT	FUNCTION
30	0	OVERALL GAIN CONTROL
	1	
	2	
	3	
	4	
	5	
	6	
	7	
31	0	NEAR GAIN CONTROL
	1	
	2	
	3	
	4	
	5	
	6	
	7	
32	0	FAR GAIN CONTROL
	1	
	2	
	3	
	4	
	5	
	6	
	7	

33	0	TD A/B	Note: Accessing this port a BMK-PULSE for the linear front end.
	1	ZBM	
	2		
	3		
	4	- INPUT PORT	
	5	-	
	6	-	
	7	-	
34	0	DRIVE SRV BY BCLK LINE: ACTIVE LOW	
	1	APPLY THE GCV SIGNAL TO LFOUT: ACTIV LOW	
	2	USE EXTERNAL GCV: ACTIV LOW	
	3	NEAR/FAR GAIN COMMUTATION: 0=EARLY 1=LATE	
	4	SRVHOLD: IF HIGH: SRV IS HELD IN THE ANALOG PART	
	5	ZINV: INVERSION BIT FOR ZBM	
	6	TRAD: CONTROLS TRANSMIT DISABLE BIT IN ANALOG PART	
	7	BLD: BEAM LOAD FOR LINEAR FRONT END	
35	0		
	1		
	2	SLOPE CONTROL FOR DYNAMIC FILTER	
	3		
	4		
	5		
	6	DYNAMIC FILTER 00: LOW 10: HIGH	
	7	FREQUENCY SELECT 01: MED 11: NO FILTER	
36	0	SHARPENING CONTROL	
	1		
	2	SELECTS CAP CIRCUIT	
	3		
	4		
	5	DYNAMIC RANGE CONTROL	
	6		
	7	EXPAND/COMPRESS	

Table 22 MEFRO.CONTROL explains the microprocessor bytes of the board.

TABLE 22

ADDRESS (HEX)	BIT	FUNCTION	
14	0	TRANSMIT FREQU. SELECT - HIGH BIT	
	1		
	2		
	3		
	4		
	5		
	6		
	7		
TRANSMIT PATTERN			
15	0	TRANSMIT FREQ. SELECT - LOW BIT	
	1		
	2		
	3		
	4		
	5		
	6		
	7		
16	0	SELECT TD A/B	
	1	RELAY 1	
	2	RELAY 2	0 = PASSIVE 1 = ACTIVE
	3	RELAY 3	
	4	TRANSMIT POWER:	0 = FULL 1 = HALF
	5	PREGAIN:	0 = FULL 1 = REDUCED
	6	BANDPASS FILTER SGM:	0 = LOW 1 = HIG BAND
	7	BOARD ENABLE:	0 = ENABLE 1 = DISABLE
17	0	CONTROL WORD FOR MICROPROCESSOR-PORT INITIALIZED WITH 80 (HEX)	
	1		
	2		
	3		
	4		
	5		
	6		
	7		

MEFRO B performs basically the same functions as MEFRO, which is used on Sigma 1SC. It is fully compatible with the old board and can be installed either on Sigma 1SC or in Sigma 1, provided that the strap ST 1 is on position REL 1.

The strap is located in the middle of the board near the bottom and the jumper must be set on the right side.

The following features are new on this board :

There is only one transmitter circuit, which is multiplexed by REL 2 either on TDA (for imaging) or on TDB (for mono-operation).


The transmitter circuit, has been modified for reduced power dissipation in the DS0026, wider bandwidth and better immunity against variations of powerfet threshold voltage.

In Sigma 1, REL 1 is used for switching the swept gain module input to the SHMD signals (TFOC output). In effect of the new ST 1 the MEFRO transmitter switches off as soon as REL 1 is activated. (If ST 1 is in REL 2 position, the transmitter can be controlled by an independent control bit).

The pattern on Z 14 and Z17 shift registers is wired with a one-bit shift compared to the old board. This delays the transmit pulse by one clock cycle (between 25 and 50 ns), but improves the timing accuracy of the transmitt signal.

3.12. Sisend p.c.b.

3.12.1. SISEND BLOCK DIAGRAM

3.12.2. DESCRIPTION OF SISEND BLOCK DIAGRAM

This board is quite similar to the ULAP SEND board. It has 16 hybrid transmitter/receiver channels and the shift/load logic required for controlling ULAP linear and curved linear transducer as shown on the SISEND Bloc Diagram. There is also a transmit-pulse generator built from a VCO and some logic followed by a programmable level translator stage in order to generate eight phase and amplitude controlled transmit control signals. The receiver outputs are connected through common-base buffers to ULAP TFOC board.

Although the functions of SISEND are quite similar to the SEND board, the microprocessor interface is quite different. It is explained in detail in a) page 58.

The hybrid circuits were also redesigned in order to put the multiplexer-control register onto the hybrid circuit. Power dissipation was also reduced. The new hybrid, which is called BIDIB, is described in b) page 60.

The voltage controlled oscillator which generates the transmitter-clock has been redesigned in order to achieve better stability. It is described in c) page 62.

There is a new transmit pattern generator which offers many possibilities in producing optimized transmit bursts.

The multiplexing of the phases for different focussing is identical to the SEND board. The power and apodization control is the same.

a) SISEND microprocessor interface

ADDR. (HEX)	FUNCTION
0	TRANSMIT FREQUENCY
1	TRANSMIT FOCUS/POWER/APOD
2	STATUS
3	TRANSDUCER CODE (INPUT)
4	-
5	-
6	ADDRESSING GENERATES MBK 
7	-

	ADDR 0		ADDR 1	ADDR 2
BIT 7 MSB	BS	BURST	POWER 2	X0 SEL LIN 1/LIN 2
6	BS 2	SELECT	POWER 1	DPLR (to TFOC)
5	BS 1	LOWER	POWER 0	BS 5 BURST SELECT
4	BS 0	BITS	APOD 2	BS 4 UPPER BITS
3	C 3		APOD 1	TEST (normally low)
2	C 2	VCO	APOD 1	FLD (fast load)
1	C 1	FREQUENCY	FOCS 1	BMR (beam rest)
BIT 0 LSB	C 0		FOCS 0	BOARD ENABLE

Oscillator-Frequency : (C0 to C3)

TRIM points F → 40 MHz
 7 → $39\frac{1}{3}$ MHz
 0 → 20 MHz

Approx. rule : $f \text{ (MHz)} = 20 + 1\frac{1}{3} \cdot C$

Burst select : (BS0 to BS5)

BS5 = low : Delay increment between transmit-phase equals VCO cycle : "lambda/8" mode.

BS5 = high : Delay increment between transmit-phase equals half VCO cycle : "lambda/4" mode.

For each case 32 possibilities.
 1F and 3F reserved for testing.

Power : (POWER 0 to POWER 2)

7 = max. power
 0 = min. power

Apodization (APOD 0 to APOD 2)

Apodization rule is : $f(n) = 1 - K \cdot \sin^2(n/16)$

APOD	K
0	2
1	1.5
2	0.92
3	0.67
4	0.45
5	0.28
6	0.13
7	0

hamming window

no apodization

Focus (FOCS 0, FOCS 1)

FOCS	TYPE
0	standard aspherical
1	far
2	close
3	curved linear

Board enable

High = enable

Shift/load logic

Clocking : address 06 HEX

Control :

FLD	BMR	FUNCTION
0	0	Load new frame*
0	1	shift 4 beam
1	0	do not use
1	1	shift 1 beam

* : works always in "wide scan" mode as shown in chapter (SCON)

Test

normal state : low

testing : high selects separate address space in the MUX EPROM.

DPLR

controls DPLR bit for TFOC as found on SEND

X0

controls optional linear transducer switchbox

X0 = low : LIN 1

X0 = high : LIN 2

Wide scanning

SISEND works always with "wide scanning".

The position of the first and the last beams relative to the transducer segments is shown in Figure SISEND WIDE SCANNING.

This way, the number of different ultrasound beams is :

$N_b = (N_c - 10) \cdot 2$ (N_c = number of transducer segments)

b) Voltage controlled oscillator

The board uses basically the same circuit as SEND. However, the control-voltage circuit has been changed. The frequency resolution has been reduced to four bits, which is enough for the application.

The basic idea of the control circuit is to directly trim three frequency points and interpolate with the DAC between them. *Figure 23 SISEND VCO* shows the basic circuit.

The most significant controlbit selects the trimming voltages through a CMOS multiplexer. The lower bits go to the DAC. For trimming the following steps are performed :

- 1 Apply the control word F (HEX). This selects the "max" trimmer for the control voltage and switches off DAC (total current on Iout).
The maximum VCO frequency of 40MHz is trimmed.
- 2 Apply the control word 7 (HEX). This selects the "mid" trimmer for the control voltage and still switches off the DAC.
The mid VCO frequency of 29 1/3 MHz is trimmed.
- 3 Apply the control word 0 (HEX). This switches the DAC on, and 7/8 of the reference current flows out of Iout.

For all other control words, the DAC will now perform a linear interpolation between the trim-points. This holds also if the MSB is high, because the "mid" trimmer is used as a DAC reference in this case.

Note that steps 1 and 2 are completely independent. Step 3 depends on 2 and should always be done last.


The circuit also has a built-in thermal compensation (by CR 1).

The stability of the trimpoints should be 3% over temperature, while the interpolated values should hold 5% of the typical values.

The control function is approximately linear.

$$f_{out} = 20 + 11/s \cdot N \text{ (MHz)} \quad \text{where } N = 0 \dots 15$$

FIGURE 23 : SISEND VCO

c) Transmit pattern generator

The transmit pattern generator is built from F 299 type shift registers. They are loaded with information coming from a PROM. The PROM can hold 64 different transmit patterns which can be different for "normal" and interlaced beams.

The pattern generator has two basically different operation modes, which are selected depending on BS5. In the 1/8 mode (see *Burst select* above) the shift registers 703/704/705 are used for pattern generation and 602 is used for transmit-phase generation. In ULAP the maximum delay is usually one wavelength. Therefore in a typical ULAP case, one bit of the transmit pattern corresponds to 1/8 wavelength.

In this mode, the delay increments of the transmit-phases equal one VCO cycle.

The delayincrement can be chosen between 25 ns and 50 ns and the "length" of a pattern-bit equals the delayincrement. Pattern length is 24 bits. In the 1/4 mode, two parallel shift register are formed by 703/704/705/601 and 701/702/703. Both are programmed with identical pattern (16 bits long). During shift operation, the first shift register chain is clocked by the true VCO-clock, and the second is clocked by the inverted clock. At the output, the taps are interlaced. This divides the effective delay increment by two. The possible delay increments range from 12.5ns to 25ns and the "length" of pattern-bit equals two delay increments.


The 1/4 mode is used for high-frequency transmit pulses only, because it offers less resolution for pattern variation.

3.13. T foc p.c.b.

The TFOC board receives the signals, originating in the transducer segments, on the eight receive phase lines.

Combinations of received phase signals are summed onto ten delay line inputs after level adjustments under program control of a ramp voltage. The summing and level adjustment provide dynamic focus of the received signal as a function of depth. Any of four delay lines can be selected by the program to match the frequency of the transducer in use.

3.13.1. TFOC BLOCK DIAGRAM


3.13.2. TFOC BLOCK DIAGRAM DESCRIPTION

The block diagram shows the main functional blocks of the board and gives a general idea of the signal flow. The computer bus extension is shown down the righthand side of the figure and it includes the eight data lines used to load data into the control registers, four address lines, and a stobe line.

a) Ramp generator

A ramp voltage, synchronous with each beam, is generated from the control word in register N° 0. The control word is applied to the DAC with its output taken to an op-amp with capacitor feedback. The capacitor is discharged through a FET gated by the send/receive signal SRV*. The DAC reference voltage is supplied by a second DAC converting the control word register N° 1 to a reference voltage. This control word provides the compensation required for curved linear transducers. The ramp generator is disabled during Doppler mode by the DPRL* signal applied to the capacitor discharge FET. A partial schematic of the ramp generator and compensation circuit is shown hereunder.


RAMP GENERATOR PARTIAL SCHEMATIC

b) Step voltage generator

This circuit uses a voltage divider followed by voltage follower op-amps in a ladder configuration to produce seven reference voltage selected by the Doppler or Imaging focus control bits of register N° 2. A partial schematic of the circuit is shown here below.


STEP VOLTAGE GENERATOR PARTIAL SCHEMATIC

c) Phase level selection circuit

This circuit, shown in the hereunder partial schematic consists of 11 comparators biased by the threshold voltage set by the Imaging focus or Doppler foccus control word. One input of each amplifier is taken from the ramp voltage and the other is connected to a step on the step voltage generator.

PHASE LEVEL SELECTION PARTIAL SCHEMATIC

d) Summing and level control

The summing and level control block shown in the block diagram uses an array of FETS with their gate level set by the output of the phase selection comparators to control the signals. The signals on the eight input lines are level adjusted and summed (the central signals of a beam are summed with a decreasing number of outer signals to give the concentration required for focus) onto ten signal lines to the inputs of the delay lines.

e) Delay line selection


The four delay lines are permanently connected to the outputs of the level control circuit and are selected via their group connections by the register N° 2 control word.

3.14. Cabling

Sigma 1 is in the great lines identical to Sigma 1SC, this is also valid for the cabling. But some points are new :

- Since there are two new boards in the analog box, there are also new cables.
(See total cabling of the two boxes in *Figure SIGMA CABLING*.)
- SISEND is connected to the ULAP-connector by a new flexible print called LIFRO. This one has to be assembled very carefully to avoid damage of the wires.
- The connection between WOMOT/MEFRO and the Wobble-Scanhead is done by INSEC. This board is redesigned, therefore its cabling has changed.

SIGMA CABLING

3.15. Video kit


See Appendix (chapter 11).

3.16. Power supply

See Appendix (chapter 11).

3.17. E.C.G. Modul

BLOCK DIAGRAM

See complete schematics in chapter 10.

ANNULAR ARRAY SUPPLEMENT

TABLE OF CONTENTS

Section 1

General information	4
---------------------------	---

Section 2

Transducers	5
-------------------	---

Section 3

Differences between Sigma 1 and Sigma 1 AC	7
--	---

Section 4

AAFRO	9
4.1 General	9
4.2 Transmission processing	10
4.3 Reception processing	11

Section 5

ATFOC	13
5.1 General	13
5.2 AA signals processing	13
5.3 AR signals processing	14
5.4 Interface	14
5.5 Ramp generator	15
5.6 Dynamic aperture comparator	15
5.7 Dynamic focusing comparators	16
5.8 Lin aperture and focus	16
5.9 AA matrix	17
5.10 LIN matrix	18
5.11 Delay lines	19


1. General information

Sigma 1 AC provides the following features :

Static aperture and focusing at the transmission


Apodization at transmission: more power is applied on the central ceramic than on the laterals to shift the secondary lobes effect.

Dynamic aperture and focusing at the reception (5 or 7 zones of dynamic aperture and 6 of dynamic focusing).

2. Transducers

The active part of AA transducers is composed of concentric 5 or 7 annular ceramics (in shape of rings).

N o t e


By convention, the ring 1 is the central ceramic and the ring 7 (or 5) is the external one.

The range of AA transducers is the following :


	AA 3.5 MHz B	
	AA 5 MHz B	
	AA 7.5 MHz B	

Another range of AA transducers is presently under evaluation and of course not released yet.

Application : A : Abdominal
B : General purpose
C : Cardiac

Sigma 1 Front end

Sigma 1 AC Front end

3. Differences between Sigma 1 and Sigma 1 AC


2 boards are completely new :

- . AAFRO (MEFRO B for Sigma 1)
- . ATFOC (TFOC for Sigma 1)
- . New cabling (MASIC)
- . New power supply
- . New keyboard
- . New AA switchbox

The equipment test (AA checkbox) was designed to check AA transducer and signals processing.


The other boards are taken from Sigma 1 :

- . SISEND II
- . SCON II
- . FIST
- . CARDIS
- . MEMEX
- . WOMOT
- . IFDOD

4.2 Transmission processing

During transmission, the MOS/FET are trimmed by TQ0-7 and TQ8-15 : one pair of FET is closed when the other one is open :


The control is obtained by the variation of the "H" stage power supply. The diodes ICR30 and 31 are protecting the first reception stage of ATFOC during transmission and REL 2 switches either to wobbler or mono transducer.

The AAFRO p.c.b. contains a VCO working identically as the VCO on SISEND II (between 20 and 40 MHz) : SRV controls the VCO running and END stops the VCO Operation.


4.3 Reception processing


During reception, IQ93 and IQ94 (triggered by the TGA) short circuit IT70. The received signal is sent to ATFOC via R3 to R9.

The HF signal (SHMD, SHMD RET) coming back from ATFOC is amplified in a sweep gain module (SGM II). The gain is adjusted according to the selected frequency and depth.

The selection of the received signal (AA, LIN or std wobbler) is done by REL 1 which switches on the appropriate channel. the LFOUT signal (HF, HFRET) is then sent to IFDOD.

ATFOC Board

5. ATFOC

5.1 General

The ATFOC board receives :


- . The eight received phase lines from SISEND (Linear transducer).
- . The seven received signals from AAFRO (AA transducer).

Each of them is sent on two separate matrixes (one for AA and the other one for linear) controlled by a zone switching controller and then delayed. The delay and summing provide dynamic focusing and aperture as a function of depth. The HF signal is sent back to AAFRO in order to be amplified.

5.2 AA signals processing


The ATFOC AA part works in two different steps :

1. Dynamic aperture.
The central ceramic is first positioned in reception (AR1) then the first and the second (AR1 + AR2) and so on to have all the ceramics in reception (AR1 to 7). This operation is carried out in the AA matrix.
2. Dynamic focusing
Dynamic focusing is obtained by applying delays on each channel. The commutation is performed by coincidence of a ramp voltage (triggered by SRV) and a comparison voltage generated by the microprocessor (according to depth, frequency and focus).

5.3 AR signals processing


Each AR signal is limited (CR) impedance matched (T) and amplified in a low capacitor differential stage which remove also the Miller effect (high Z_i , low Z_o) before being sent on the AA matrix.

The 5/7 signal earthes the AR1 and 7 when operating with a 5 rings AA transducer.


5.4 Interface

The data bus (d0 to d7) is at first sent to a buffer (Z8) and then to four A/D converter (Z1 to 4) in order to generate the comparison levels and command signals used to trigger the matrix. Each A/D converter has four outputs which are adapted by an amplifier stage (Z1 to Z8). L1 to L13 are the comparison voltages selected according to the chosen depth. Z10 is a latch which is used as a pointer for the A/D converter. Then Z10 selects by DLS1 to 3 the delay line(s) according to the transmitting frequency.


5.5 Ramp generator

The ramp generator is performed by an Integrator (Z13) triggered by L14. The integrator circuit is composed of Z13, C3 and R119 and the output signal is limited to 5.5 V by Q74 and R118. SRV signal discharges C3 at the transmission end to begin the ramp generation. RAMPS and RAMPL are sent to the sector and linear matrixes (TP16 and 17).


5.6 Dynamic aperture comparator

The L1 to L7 signals are compared to RAMPS in the comparator stages (Q22 to Q42). When one of the input Ln is in accordance with the programmed level, the corresponding transistor let the SAn signal flow to the AA matrix to open the corresponding ceramic.


5.7 Dynamic focusing comparators

As soon as the last ceramic is opened (L7, SA7), Q49 short circuits the dynamic aperture comparators (except SA4 which is always active).


By connecting together the collector of the different stage transistors (ex : Q8 and Q7), SB2 is activated by the \bar{V} of SB1, etc. ...

5.8 Lin aperture and focus

First L1 and L2 are sent to the comparator stage (Q68, Q72) to generate the aperture orders (LA1, LA2). Then focusing is applied (L3 to L8, Q50 to 66, LB1 to LB6).

5.9 AA matrix

The AA signals are first adapted and limited. 5/7 (L16) signal earthes AR1 and AR7 in case of using a 5 rings ceramics. SA1 to 7 activate the Q38 to 43 transistors for the aperture. Then SB1 to 6 selects the received channels on which delay has to be applied by the delay lines (N0 to N18) for the dynamic focusing. In order to have a smooth switching commutation inside the matrix, the AR signals are divided by 3 (R83 to R104).


		Element No (SAR1 to SAR7)						
Focus range		1	2	3	4	5	6	7
SA1	1	18	-	-	-	-	-	-
SA2	2	18	15	-	-	-	-	-
SA3	3	18	15	12	-	-	-	-
SA4	4	18	15	12	9	-	-	-
SA5	5	18	15	12	9	6	-	-
SA6	6	18	15	12	9	6	3	-
SA7	7	18	15	12	9	6	3	0
SB1	8	15	13	11	9	7	5	3
SB2	9	12	11	10	9	8	7	6
SB3	10	9	9	9	9	9	9	9
SB4	11	6	7	8	9	10	11	12
SB5	12	3	5	7	9	11	13	15
SB6	13	0	3	6	9	12	15	18

Element	1	2	3	4	5	6	7	8	9	10	11	12	13
1													
2													
3													
4													
5													
6													
7													

5.10 LIN matrix

LB1 to LB6 activate the LIN Matrix commutation. PRO to 7 (from SIDEND) are sent to the delay lines (DL1 = 199.8 ns, DL2 = 333 ns, DL4 = 468.9 ns). In order to have smooth switching commutation inside the LIN matrix, the RP2 to 7 signals are divided by 3 (R196 to R111). RP0 and 1 are used for dynamic aperture (2 zones).

		Element No							
Focus		1	2	3	4	5	6	7	8
25 mm	1	0	30	120	—	—	—	—	—
32 mm	2	0	30	90	180	270	—	—	—
50 mm	3	0	30	60	120	180	270	390	—
70 mm	4	0	0	30	90	120	210	270	360
100 mm	5	0	0	0	30	60	90	120	180

5.11 Delay lines


DLS1 and DLS3 signals enable the delay lines according to the dynamic focusing .
 N0 to N8 at the delay line inputs are the outputs of the AA matrix (if RAMPS =1) or LIN matrix (if RAMPL=1).
 LIN signal short-circuits N3,N5,N15 and N18 in case of linear transducer, and SEC short-circuits N4 and N1 in case of 5 rings AA transducer.

Note that N0 is the input without delay to be applied.

N0 to N18 are delayed and summed in DL1 to DL3 (TP15), then impedance adapted (T8) and finally gain adapted (TP 16). SHMD is the output signal to be send to AAFRO.


NOTE

To operate in right conditions, the delay line must be terminated by a 100 Ohm load.


The Hi inputs are the common stages Q104 to Q120. The 100 Ohms input loads are R132, R134 and R136. The output 100 Ohms loads are R131, R133 and R135 (205 Ohm) in parallel with T8 (205 Ohm).


The outputs of the delay line are voltage limited before triggering the adaptation stage.

TRANSUDERS CODDING

ANNULAR TRANSDUCER

WOBBLER TRANSDUCER

PART NUMBER	TRANSDUCER TYPE	STRAP ON
597 430	Annular Array Transducer 3.5 MHz type A	PR1 PR2 PR4
594 970	Annular Array Transducer 3.5 MHz type B	PR0 PR2 PR4
597 740	Annular Array Transducer 5 MHz type A	PR0 PR1 PR2 PR3
592 080	Annular Array Transducer 5 MHz type B	PR0 PR1 PR2 PR3
597 120	Annular Array Transducer 7.5 MHz type B	PR3 PR4

PART NUMBER	TRANSDUCER TYPE	STRAP ON
856 150	Wobbler Transducer 3.5 MHz type A	PR0
581 410	Wobbler Transducer 3.5 MHz type C	PR1
589 330	Wobbler Transducer 3.5 MHz type D	PR0 PR1 PR3
856 010	Wobbler Transducer 5 MHz type B	PR2
583 170	Wobbler Transducer 7.5 MHz type B	PR0 PR2
598 880	Wobbler Transducer 7.5 MHz type V	PR0 PR3 PR4

Adjustment and Calibration

4. ADJUSTMENT AND CALIBRATION (UNIT OPEN)

4.1. Introduction

This chapter describes the adjustment and calibration of an open unit. It is mandatory to perform a complete unit check and a safety check afterwards (See chapter 5).

Functional safety is ensured by performing the relevant test in a professional manner.

If the unit is opened, it is imperative to perform all safety and functional checks (described in chapter 5 of this manual) after reassembly, to ensure that functional safety demands are met.

Caution : CMOS components : See chapter 2.2. page 5.

4.2. Calibration equipment

4.2.1. TEST OF RECEIVE ELECTRONIC (common to linear and sector)

*You need : Open machine
Sigma I Sector Checkbox
Oscilloscope*

- Connect the transducer cable as well as the TM/Doppler connector of the checkbox to the Sigma 1.
- Select the following settings on the checkbox :
 - ECG off
 - XMIT off
 - MODE 4
 - PULSE LENGTH F
 - TRANSDUCER CODE FD (= Wobbler 3.5MHz C)
 - FREQUENCY 3.5MHz

a) Control of TGC


- Press TM
DEPTH 4
FAR
- Set TGC1 ... TGC 9 at maximum
- Observe the GCV voltage (available on the checkbox) on the oscilloscope (1V/div, 50 micro-sec/div, synchronized on SRV  (available on the checkbox)).
- Control on the scope that $GCV = 4.0V \pm 200mV$ for TM gain min
 $= 5.0V \pm 200mV$ for TM gain max
- Verify that the time which elapses between  and the beginning of GCV reset is included in the interval 239...248 us.
- Set TCG 2, TCG 3 ... TCG 9 at minimum.
- Adjust TCG 1 to suppress the break of the first zone.
- GCV must be absolutely constant. Any slope must be corrected on the CARDIS board with R 166.
- Control on the scope that $GCV = 0.60V \pm 0.1V$ for TM gain min.
 $= 3.60V \pm 0.1V$ for TM gain max.

b) Control of the log amp tuning and overall gain calibration


- All other settings being maintained, set :
TM max
Reject max
- Observe the TMZ output (available on the checkbox) on the scope (500mV/div, 50 micro-sec/div synchronized on SRV—).
According to the jumper on the CARDIS board, there are 2 possibilities :

	POSITION B	POSITION A
White level	$2.77V \pm 100mV$	$3.40V \pm 100mV$
Black level	$1.02V \pm 100mV$	$0.00V \pm 100mV$
Blank level	$0.40V \pm 100mV$	$-0.25V \pm 100mV$

- Set Reject min
Enhance 0
- Adjust the TM Gain to have the max of TMZ at 2.55V (2.84 V with jumper in position A).
- Observe the GCV voltage on the scope (500mV/DIV, 50 micro-sec/div, synchronized on SRV —). It should be at $2.00V \pm 200mV$.
- Select the mode 0 of the checkbox and observe again the TMZ output on the scope (500mV/div, 50 micro-sec/div, synchronized on SRV—).
You should see 9 equally spaced curves and the top of the 10th curve should be almost distinguishable.

Jumper in position B


Jumper in position A

- Select again the mode 4 of the checkbox.
- Choose Enhance -1 (= ACG on).
- The top of the curve should be at $2.0V \pm 100mV$ (2.12V with jumper in position A).

c) Control of Tracking Filter Tuning


- All the setting being maintained, observe the image on the SIGMA 1 TV monitor. The center of the max must lie at about 9.5cm (± 0.5 cm).

- Select Frequency 5MHz on the checkbox.
- Press DEPTH 1
NEAR
- Observe the TV monitor. The max must lie at about 6.4 cm.


d) Control of the Transfer of the Digital Part

- Set normal TM (white on black)
- Press TM
POST 4
SET T GREYSCALE ALL
- Observe the TV monitor. Make sure that all the bits of the AD converter are present (smooth picture).
- Control the dynamic range of the AD converter on the FIST board. With a digital voltmeter you must find 3.00V (± 0.01 V) at TP 10.


e) Control of the pulse Shape

- Select the following settings on the checkbox :
 MODE 6
 PULSE LENGTH F
 FREQUENCY 5MHz
 (ECG off)
 (TRANSDUCER CODE FD)
- Press TM
 DEPTH 4
 NEAR
 SET T DIRECT-TMZ
- Press Enhance = 0
- Adjust TM gain at 1.7V (scope at 500 mV/div, 50 micro-sec/div, synchronized on SRV ).
- Observe the TV monitor and place the pulse at depth of 5cm with the DEPTH potentiometer of the checkbox.
- Observe the TMZ voltage on the scope (500 mV/div, 50 micro-sec/div).
 The 8 curves should look like those of the photo below, without over-or undershot and with fast rise and fall time.

Jumper in position B

Jumper in position A

4.2.2. CONTROL OF THE SEND FREQUENCY TUNING


*You need : Open or closed machine
Sigma 1 Sector Checkbox
Oscilloscope
50 Ohm termination*

a) Control of the frequency tuning


- Select XMIT ON on the checkbox.

Connect the transmission output of the checkbox to an oscilloscope with a coax cable and 50 Ohm termination. Observe the pulse on the scope (500mV/div, 100 ns/div, synchronized on SRV \rightarrow). It must be clean without overshoot, undershoot or oscillating tendency. Within, the measurement accuracy you should find the indicating time interval between the two markers.


- Press 2D
DEPTH 3
NEAR

- Press DEPTH 3
FAR


- Press DEPTH 4
NEAR


b) Control of the pulse amplitude

- Press DEPTH 3
NEAR

- Measure the peak to peak voltage of the pulse. Within $\pm 0,1V$, you should find the indicated one.

- Press DEPTH 1

c) Control of the monochannel access

- Connect the monochannel output of the Sigma 1 to an oscilloscope with a coax cable and 50 Ohm termination.
- Press DEPTH 3
NEAR
TM
SET TM MONO
- Observe the pulse on the scope (20V/div, 200 ns/divsynchronized on SRV_). It should look like the one depicted on the previous page.


4.2.3. TEST OF LINEAR TRANSMIT AND FOCUSING ELECTRONIC

*You need : Open or closed machine
Sigma 1 linear checkbox
Oscilloscope*

- Connect the transducer cable as well as the TM/Doppler connector of the checkbox to the Sigma 1.
- Select the following settings on the checkbox :
XMIT ON
MODE 9
TRANSDUCER CODE 7A (= ULAP lin 3.5MHz)
FREQUENCY 3.5MHz
PULSE LENGTH F
SEGM 6


a) Control of the frequency tuning

- Connect the transmission output of the check box to an oscilloscope with a coax cable and 50 Ohm termination. Observe the pulse on the scope (500mV/div, 200 ns/div, synchronized on SRV \rightarrow). Within the measurement accuracy you should find the indicated time interval between the two markers.
- Press TM
DEPTH 5
FAR

b) Control of the transmit focussing

All settings being maintained, measure on the scope the pulse amplitude. You should find $0.95V \pm 0.1V$. By selecting the different segments on the checkbox, make sure that all 16 phases are equal (= within the tolerance above specified) and with the same polarity. If you find an idle phase, you have got an interlace beam. In this case, move the TM marker 1 step with the trackball and make sue that all 16 phases are working.

- Press NEAR


By selecting the different segments, control the transmission apolization.

c) Control of the receive-focussing

All settings being maintained,

- press 2D
DEPTH 4
NEAR

Observe the TV monitor. The picture must look like the photo below. By selecting the different segments on the checkbox, make sure that all the 16 receive channels are OK.

d) Control of the linear Gain Calibration

- Select the mode 0 of the checkbox
- Press 2D
DEPTH 4
NEAR
- Place the TM marker on the brightest bar
- Press TM
- Set TGC2, TGC3, ... TGC9 min
Repeat min
ENHANCE 0
- Observe the GCV voltage (available on the checkbox) on the scope (500mV/div, 50 micro sec/div, synchronized on SRV \rightarrow).
- Adjust TGC1 to suppress the break of the first zone.
- Observe the TMZ output (available on the checkbox) on the scope (500mV/div, 50 micro sec/div, synchronized on SRV \rightarrow).
- Adjust the TM gain to have the max of TMZ at 2.00V (1.90V with jumper in position A).
- Observe again the GCV voltage on the scope. It should lie at 1.80V ± 200 mV).


e) Control of the pulse slope

All other settings beeing maintained,


- Select the mode 6 of the checkbox
- Press 2D
- Place the TM marker on the brightest bar
- Press TM

SET T DIRECT-TMZ

- Observe the TV monitor and place the pulse at a depth of 9.5cm with the depth potentiometer of the checkbox.
- Observe the TMZ voltage on the scope (500mV/div, 50 micro sec/div, synchronized on SRV $\overline{\text{f}}$, delayed trigger 114 micro sec and 2 micro sec/div). The curve should look like the one on the photo below, without over or undershoot and with fast rise and fall time.

Jumper in position B

Jumper in position A


4.2.4. TEST OF THE VIDEO OUTPUT

*You need : Open machine
3.5MHz cardio probe
75 Ohm termination
Oscilloscope*


- Connect a 3.5MHz cardio transducer
- Set normal image
 - Post 2
 - TGC 1, TGC 2 ... TGC 9 min
 - 2D Gain min
 - TM Gain min
 - Reject min
- Observe the video signal of the Rec Out output terminated with a 75 Ohm resistor on the scope (200mV/div, 10 micro sec/div).
- With the black porch as reference level, you should find the following voltage measurements within $\pm 50\text{mV}$:

Overlay extra-white	710mV
Image white for normal image	700mV
Image white for inverse image	480mV
Overlay white	480mV
Image black	40mV
Black porch	0mV
Sync pulse	-300mV

If there is an oscillation, the 16 shades of the greyscale are not clearly visible and the image black is thick. This can better be seen in inverse TM (black on white). The two photos show a normal Rec Out output.

2D

TM

- Control of the NSTC Standard :
Switch power off. Set the jumper ST2 of the PERIP Board on NTSC and the TV monitor switched on 60Hz. The geometry must still be OK. Return to the PAL standard (PERIP ST2 on PAL and TVM switch on 50Hz).

4.2.5. TEST OF THE LSR OUTPUT

*You need : Open machine
3.5 cardio probe
Line scan recorder*

- Connect a 3.5MHz cardio probe and make an image of the RMI phantom.
- Connect the ECG simulator.
- Press ECG
TM
SET TM 50mm/s
SET T 16 GREYSCALE
DEPTH 1
- Record on paper about 5 seconds for each depth.
After each DEPTH-change you must press again on the sofkey 16 GREYSCALE.
The time-depth markers must look like black dots surrounded with a white rectangle. Verify that the depth indicated by these markers correspond to the nominal one within 4%.
Check the quality of the ECG trace.
- Press 2D
- Place th TM beam on a small detail of the RMI phantom.
- Press TM
SET TM 25 mm/s
- Record on paper about 5 seconds for each depth and check the quality of the record.
- Control of the PRINT key :
Connect the REC OUT output of Sigma 1 to the video input of the line scan recorder and make a copy using the PRINT key of Sigma 1.

4.2.6. TEST OF DOPPLER HF INPUT/OUTPUT SOCKET


The connection between the Doppler HF socket on the rear panel and the connectral have to be ckecked with a special test adapter (to be designed) to make sure that the switching on the INSEC board works properly and that there is no crossing between the signals (the crossing of Tx+ and Tx- for instance has very undesirable effect).

4.2.7. CONTROL OF THE INTERFERENCE LEVEL

*You need : Closed machine
3.5MHz cardio probe
5MHz ULAP linear transducer*


- Connect a 3.5MHz cardio probe and a 5MHz ULAP linear transducer to the Sigma 1 closed.
- Set TGC 1, TGC 2, ... TGC 9 max
2D Gain max
TM Gain max
Reject min
- Pres TM Normal (white on black)
- Set Trad = 0 (with the switch on the Doppler/TM connector)
- Observe the Sigma 1 TV monitor. The noise appearance must be steady and smooth, without any distinguishable dot, line, curve or pattern. Try it with the 5 depths.
- Press 2D Normal (white on black)
Idem
- Press LIN
Idem

4.3. TV monitor p.c.b.


TP 21

0V
H
0,5V/DIV


TP 5

0V
H
1V/DIV


TP7

H
10V/DIV
0V




TP4 (1)

H
0,5V/DIV
0V

TP4(2)

V
0,5V/DIV
0V

Maintenance and Safety Checks

5

5. SAFETY AND PERFORMANCE CHECKS (UNIT CLOSED)

5.1. Test of closed machine

KEY	ACTION
SET TECH DATA	Softkey-line changes
SOUND VEL	Sound velocity 1540 appears. Program 1401
SET CLOCK	Program permanent clock
ENHANCE	Return to standard softkey-line
FREEZE	Defreeze
DEPTH 1	Image goes to 8 cm depth
ANGLE 60°	Angle changes to 60°
SET TECH DATA	Softkey-line changes
INIT	Stores init-format into permanent RAM
	Switch power off and wait 10 seconds. Switch power on Check if depth is 7.5 cm, angle 60°, date as programmed above.
SET TECH DATA SOUND VEL	Check if sound vel. is 1401. Program 1540
	Note : Loss of angle, depth, date or sound velocity indicates a problem with the permanent RAM. Check battery.

KEY	ACTION
FREEZE	Defreeze
ENHANCE	Return to standard softkey-line

Now different image formats should be tested. There are 48 different formats. We do not test all of them but just the extreme ones.

Check if machine is in 60°, 8cm, high resolution. For the format above, look if the image is stable, without white or black lines. Step through the following states:

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KEY	ACTION
SET I	Brings you back to 60°, 8cm
NEAR/FAR	Switch back and forth from near to far. Near mode should give a low gain image with very fine granularity. Far mode should give a bright image.
FREEZE	freeze image
FREEZE	image goes back to live
Footwitch FREEZE	live image
→ ←	inverts image orientation
TM	format changes to TM
SET TM	softkey-line changes
10 50 20	check different TM speeds
MONO	Blank TM trace if no MONO Transducer connected. Connect one and press it on the phantom. Check result.
TM	return to wobbler
SET	fetch test menu
T	
GREYSCALE ALL	generates all 60 greylevels. Checks homogeneity. No steps may be visible.

KEY	ACTION
Footswitch 2D/TM	TM-image
Footswitch 2D/TM	2D-image
ENH	try enhance -1, -2, -3 characteristic changes are visible
ENH	step through -2, -3, 0, 1, 2, 3 characteristic changes are visible. Go back to ENH 0.
PROC	step through PROC 2, 3, 4 characteristic changes are visible
FILTER	step through FILTER 1, 2, 3 small changes in signal to noise are visible. Compare with reference machine.
SET I	go back to init states connect ECG-simulator to ECG input. Switch on standard ECG.
ECG YES/NO	ECG appears on screen, must be clean switch all interference sources of ECG- simulator on (50 Hz, base line shift) ECG on screen will be less clean but still well visible
GAIN+	Check ECG gain control.
GAIN-	Adjust to full scale signal
SYNC I	check sync mode use soft keys for adjusting sync positions
TM	TM-image
ECG YES/NO	switch ECG off
2D	2D-image use trackball for adjusting TM-marker

KEY	ACTION
TEXT	use trackball for adjusting cursor
A-Z	set cursor up - left. Try all letters Q...M and 1-9
1-9	numbers 1...9 and SPACE - / . LINE FEED
	check rub-out
CLEAR	clear screen (takes several seconds)
SET	cursor jumps to PAT ID-field
PAT ID	write QUERTY
PAT ID	clears PAT ID
PAT ID	recalls "QUERTY"
SET	
TECH DATA	SIGMA 1 write permanent TEXT : KONTRON into the upper left corner exit with "LINE FEED" "LINE FEED"
TECH DATA	recall, clear, recall again the permanent text
TECH DATA	clear and recall technical data
MAG	magnifies the image; use trackball for positioning the window
MAG	return to normal mode
DIST	check distance measurement on graticule
\pm	
AREA	check area measurement
\pm	
PRINT	makes a hard copy on the line scan recorder (if there is any, else measure output on rear panel connector) Now you have controlled all keys and all operating modes of the machine. Now the correct time and date and the init mode should be programmed.

KEY	ACTION
SET I	recalls actual init mode
80° DEPTH 3	corrects for wanted mode
SET TECH DATA MORE INIT	stores init mode

END OF TEST

Test of external TV connections

Conect TV monitor with 75 Ohm impedance to "Monitor out".
The image should appear as on the monitor.

Switch live-replay to REPLAY. Both images should disappear.

Connect TV-monitor with 1k impedance to "REC out" and to "REC in".
Now the image appears again on both monitors.

Switch to LIVE.

Now the intensity of the external monitor must **not** change. If it changes, there is something wrong with the dummy 75 Ohm termination on the live-replay switch.


PS : Make sure that every machine leaving the factory has the following
init-mode setting :

- Depth 4
- HR (high resolution)
- NEAR (near mode)
- Angle 80°
- Enhance off
- Filter off
- Image polarity : normal (black background)
- TM polarity : inverse (white backgroud)
- Sound velocity : 1540 m/s
- Correct date
- Processing 3
- Permanent Memory enmpty

5.2. Safety tests

Note : the measuring frequency of the ampere meter may be reduced to 50kHz with a low pass filter.

TEST I

D.U.T. : device under test

ECG : ECG input

P : phase

N : neutral

GND : protection ground

I_L : leakage current


Condition : $I_L < 100 \text{ micro-A}$

TEST II

Same as I but neutral disconnected

Condition : $I_L < 100 \text{ micro-A}$


TEST III

220V AC applied to any patient ECG lead

Condition $I_L < 50$ micro-A

TEST IV


GND disconnected

Neutral connected or hot (whichever is worse)

Condition : $I_L < 50$ micro-A

TEST V

Normal operation

Condition : $I_L < 10 \text{ micro-A}$

Assembly and Disassembly

6

6. ASSEMBLY AND DISASSEMBLY

6.1. CMOS components

The SIGMA 1 has CMOS parts. During maintenance it is essential to prevent transfer from static electrical charges to circuit board and components by :


- grounding instrument
- grounding soldering iron
- wearing a grounded wrist strap (via 1 Mohm to ground)

It is imperative to keep the surface of the boards absolutely clean and insulated.

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








Trouble Shooting

7

7. FAULT FINDING AND STANDARD EXCHANGE

* in preparation *

Parts List

8. PARTS LIST

8.1. Recommended spare parts list

PART NUMBER	DESCRIPTION	A	M.B	B	C
	BOARDS (incl. components)				
	=====				
852724	IFD00	*	1	2	0
852775	MEFRO + SGM2	*	1	2	0
854956	PERIP + PARSIO + FLAT CABLE	*	0	1	0
861286	SCON 2	*	1	2	0
852848	MEMEX	*	1	2	0
852864	FIST	*	1	2	0
852899	WOMOT	*	1	2	0
852929	CARDIS	*	1	1	0
854786	COMPHO	*	0	1	0
863491	INSEC B	*	0	1	0
854751	POLIN	*	0	1	0
861308	SISEND	*	1	2	0
845604	T-FOC	*	1	1	0
852791	LIFRO	*	1	3	0
	SUBASSEMBLIES				
	=====				
861561	TV MONITOR (PCB & CRT)	*	0	1	0
854638	POWER SUPPLY	*	1	1	0
855529	TRACKBALL	*	1	1	0
855413	ECG AMPLIFIER	*	0	1	0
862827	KEYBOARD SIGMA 1	*	0	1	0
860859	CLAVIER FRANCAIS	*	0	1	0
	CABLES				
	=====				
858781	FLAT CABLE N 4	*	0	1	0
858803	FLAT CABLE N 6	*	0	1	0
858811	FLAT CABLE N 10	*	0	1	0
858838	FLAT CABLE N 11	*	0	1	0
858862	FLAT CABLE N 18A	*	0	1	0
858897	CABLE FOR TRACKBALL N18D	*	0	1	0
856894	ECG CABLE	*	0	1	0
855049	WOBBLER TRANSDUCER CABLE	*	0	1	0
	DIVERS				
	=====				
	BOTTOM COVER	*	0	1	0
857998	FRONT COVER (KEYBOARD HOUSING)	*	0	1	0
860131	KONTRON INSTRUMENTS LABEL FOR FRONT COVER	*	1	3	0
856339	FRONT PANNEL	*	0	1	0
860832	FACE DECOR	*	0	1	0
857106	COVER	*	0	1	0
860611	SECTORIAL TRANSDUCER HOLDER	*	0	1	0
847275	LINEAR TRANSDUCER HOLDER	*	0	1	0
850608	ROLL FOR CAMERA	*	2	3	1
	TEST TOOLS				
	=====				
861448	TEST CABLE	*	1	1	0
587370	SECTORIAL CHECK BOX	*	1	1	0
586730	LINEAR CHECK BOX	*	1	1	0
	DOCUMENTATION				
	=====				
867381	SIGMA 1 SERVICE MANUAL	*	1	2	1
864056	OPERATING MANUAL	*	1	1	1
864064	MANUEL D'UTILISATION	*	1	1	1

[illegible]

RECOMMENDED SPARE PARTS LIST FOR SIGMA 1 AC

PAGE 1 OF 2

PART NUMBER	DESCRIPTION	A	M.B	B	C
	BOARDS (incl. components)				
	=====				
868116	AAFO	*	1	2	0
854956	PERIP + PARSIO	*	0	1	0
861286	SCON 2	*	1	2	0
852848	MEMEX	*	1	2	0
852864	FIST	*	1	2	0
852899	WOMOT	*	1	2	0
852929	CARDIS	*	1	1	0
854786	COMPHO	*	0	1	0
868167	MASIC	*	0	1	0
854751	POLIN	*	0	1	0
861308	SISEND	*	1	2	0
868132	AT-FOC	*	1	1	0
852791	LIFRO	*	1	3	0
	SUBASSEMBLIES				
	=====				
861561	TV MONITOR (PCB & CRT)	*	0	1	0
874914	POWER SUPPLY AC	*	1	1	0
855529	TRACKBALL	*	1	1	0
855413	ECG AMPLIFIER	*	0	1	0
874523	KEYBOARD SIGMA 1 AC	*	0	1	0
860859	CLAVIER FRANCAIS	*	0	1	0
	CABLES				
	=====				
858781	FLAT CABLE N 4	*	0	1	0
858803	FLAT CABLE N 6	*	0	1	0
858811	FLAT CABLE N 10	*	0	1	0
858838	FLAT CABLE N 11	*	0	1	0
858862	FLAT CABLE N 18A	*	0	1	0
858897	CABLE FOR TRACKBALL N°18D	*	0	1	0
856894	ECG CABLE	*	0	1	0
	MISCELLANEOUS				
	=====				
	BOTTOM COVER	*	0	1	0
857998	FRONT COVER (KEYBOARD HOUSING)	*	0	1	0
860131	KONTRON INSTRUMENTS LABEL FOR FRONT COVER	*	1	3	0
856339	FRONT PANNEL	*	0	1	0
860832	FACE DECOR	*	0	1	0
857106	COVER	*	0	1	0
860611	SECTORIAL TRANSDUCER HOLDER	*	0	1	0
847275	LINEAR TRANSDUCER HOLDER	*	0	1	0
850608	ROLL FOR CAMERA	*	2	3	1
	TEST TOOLS				
	=====				
861448	TEST CABLE	*	1	1	0
587370	AA CHECK BOX	*	1	1	0
586730	LINEAR CHECK BOX	*	1	1	0
	DOCUMENTATION				
	=====				
867381	SIGMA 1 SERVICE MANUAL	*	1	2	1
594310	OPERATING MANUAL	*	1	1	1

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005

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Edi. Da

Product Group : Ultrasound
Ref. N° : 87190 K.I.S. BM/mm
Date : October 14th, 1987
Contact : Bernard MARTIN

Sigma line Engineering Service Information N° 1

**Subject : new soft compatibility level 6 for Sigma 1 SC
Sigma 1
Sigma Doppler**

Four new softs are available :

- . 4 proms V0.4/C6 to be implemented on SCON 1 for Sigma 1 SC
- . 5 proms V1.11/C6 to be implemented on SCON 2 for Sigma 1 LIN
 Sigma 1 STAR
 Sigma 1 CLASS
- . 5 proms V3.11/C6 to be implemented on SCON 2 for Sigma 1 CARDIO
- . 5 proms V0.9/C6 to be implemented on SCON 2 for Sigma Doppler

Caution

**If Sigma Doppler is part of the system to be
upgraded, both Doppler soft and Sigma soft must
be exchanged in order to ensure compatibility.**

**Be sure you use adequate chip type when
duplicating the proms.**

Contents of this Tech Bulletin :

Page 2	How to check the soft version of the unit
Page 3	Proms exchange
Page 4	Prom types and part numbers
Page 5	Main improvements after implementing soft C6 on Sigma 1 SC
Page 6	Main improvements after implementing soft C6 on Sigma 1
Page 7	Main improvements after implementing soft C6 on Sigma Doppler

How to check the soft version implemented in the unit

Sigma 1 SC

Press the [SET] [V] touchkeys and the actual soft version appears in the right upper corner on the CRT. After implementing the new soft, the display will show :

V0.4/C6

Sigma 1

Press the [SET] [V] touchkeys. After implementing the new soft, the display will show :

V1.11/C6

for Sigma 1 CLASS
Sigma 1 STAR
Sigma 1 LIN

V3.11/C6



for Sigma 1 CARDIO

Sigma Doppler


Press the [SET] [MORE] [VERSION] touchkeys. After implementing the new soft, the display will show :

V0.9/C6

Prom exchange

CAUTION : Strap inverted →

Prom types and part numbers

Unit	Prom P.N.	Specification	Location on SCON	Type description
Sigma 1 SC	858 668	630 032	Z 110	27128 K
	858 641	630 033	Z 111	27128 K
	858 633	630 034	Z 611	27128 K
	858 617	630 036	Z 614	2764-30
Sigma 1	866 903	630 082	Z 110	27128 K
	866 881	630 081	Z 109	27128 K
	866 911	630 083	Z 512	27128 K
	866 938	630 084	Z 513	27256
	866 946	630 085	Z 514	27256
Sigma Doppler	867 764	630 094	Z 110	27128 K
	867 772	630 093	Z 109	27128 K
	867 799	630 095	Z 512	27128 K
	867 802	630 096	Z 513	27128 K
	867 829	630 097	Z 514	27256 or 27C256

Main improvements after implementing soft version V0.4/C6 on Sigma 1 SC

- . New TM speeds : 25 mm/s - 50 mm/s - 100 mm/s (Instead of 10 - 25 - 50)
- . Depth markers will be displayed in TM.
- . Scanning depth 1 using 7.5 type B transducer will be 6 cm (Instead of 7.5 cm).
- . Automatic freeze delay is protracted from 5 minutes to 15 minutes.

Main improvements after implementing soft version V1.11/C6 or V3.11/C6 on Sigma 1

- . Direct access to the biometric functions by pressing [SET] [C] touchkeys.
- . New TM speeds : 25 mm/s - 50 mm/s - 100 mm/s (Instead of 10 - 25 - 50)
- . Simpson and hemi-ellipse : the systole volume and diastole volume can be computed.
- . The technical data are automatically displayed when exiting the CINE mode.
- . Set up programs : the filters automatically selected are changed :

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Main improvements after implementing soft version V0.9/C6 on Sigma Doppler

- VECTOR : Maximum vector angle is 70° (Instead of 45° max).
- TM : An additional speed of 100 mm/s is added.
- BIOMETRY : New TM measurements are implemented :
DIST
RATIO
SLOPE
TEICHHOLZ
- BIOHELP : During a Doppler measurement function, a text used to help the operator is displayed on the screen : [SET] + [MORE] + [BIO-HELP] + [ON] .
- CLEAR : Enables the user to erase a Doppler biometry curve.
- GATE : Between two measurements, the Doppler gate position can be modified by using the [TRACKBALL] by pressing [DEFREEZE]. The gate position can also be modified in Magnifier mode by using the [TRACKBALL] after pressing the [POS TM] touchkey on Sigma 1.
- 2D : measurements on 2D picture are possible.
- 2D and TM REFRESH : The 2D or TM picture on the Sigma 1 screen can be refreshed during Doppler examination : [SET] function ([SET] + [2D REFR] + [ON]).

Product Group : Ultrasounds
Ref. N° : 87203 K.I.S. BM/mm
Date : October 14th, 1987
Contact : Bernard MARTIN

Sigma 1 & Sigma Doppler
Engineering Service Information N° 2

Problem : Unit self resets while operating without any apparent reason, on units equipped with the C6 level soft version.

Modification : Cut pin 11 of Z511 (LS373). Connect this pin (using a short wire) with pin 13 of Z405 (LS02).

Priority code : 3
The modification may apply only if the problem actually appears on installed units. According to our experience, this problem should occur very rarely, assumed that the C6 soft version is implemented.

Attn.: _____

Product Group : Imaging
Ref. N° : 87259 K.I.S. BM/mm
Date : December 15th, 1987
Contact : Bernard MARTIN

Sigma 1

Engineering Service Information N° 3


Problem	Unusual noise in the sector image or in the linear image (especially if using 5 MHz linear transducer).
----------------	---

Cause	The surface of ANREF shielding (between logic box and analog box) has lightly become corroded, and therefore avoids perfect grounding.
--------------	--


Solution	See page 2 of 2
-----------------	-----------------

Modification	Priority code is : 3 This modification must be performed on units when described problem occurs, or if a unit is at hand for repair or during preventive maintenance.
---------------------	--

Ordering information	Requested part (folded spade terminal - ref.: 842 613) is available f.o.c. on request to myself.
-----------------------------	--

- 1° Remove the ANREF p.c.a.
- 2° Clamp the two contact parts.
- 3° Clean ANREF and both sides of aluminum with alcohol.
- 4° Remount ANREF and screw the two nuts.

Product Group : Imaging
Ref. N° : Sigma Line - E.S.I. N° 4
Date : April 18th, 1988
Contact : Bernard MARTIN

Sigma Line

Engineering Service Information N° 4

Subject : updated software for Sigma 1, Sigma 1 AC

Seven new softs are available :

. 5 proms	V1.12/C6	for Sigma 1 LIN for Sigma 1 STAR for Sigma 1 CLASS
. 5 proms	V3.12/C6	for Sigma 1 CARDIO
. 5 proms	VS1.12/C6	for Sigma 1 AC STAR for Sigma 1 AC CLASS
. 5 proms	VS3.12/C6	for Sigma 1 AC CARDIO
. 5 proms	VA1.1/C6	for Sigma 1 AA STAR
. 5 proms	VA3.1/C6	for Sigma 1 AA CARDIO
. 5 proms	V0.11/C6	for Sigma DOPPLER

Caution

This Technical News Letter does not concern Sigma 1 SC. There are no changes regarding its compatibility level (C6).

Make sure that you use correct chip type when duplicating the proms.

One set of proms will be despatched f.o.c. to each country on written request (telex to B. Martin), mentioning the desired software version(s).

Contents of this Technical News Letter :

Page 2	How to check the software version of the unit
Page 3	Proms exchange
Page 4	Prom types and part numbers
Page 5	Main improvements after implementing software C6 on Sigma 1, on Sigma 1 AC, on Sigma 1 AA, and on Sigma Doppler

How to check the soft version implemented in the unit

Sigma 1

Press the **SET** **V** touchkeys and the actual software version appears in the right corner on the CRT. After implementing the new software, the display will show :


V1.12/C6	for	Sigma 1 CLASS Sigma 1 STAR Sigma 1 LIN
V3.12/C6	for	Sigma 1 CARDIO
VS1.12/C6	for	Sigma 1 AC STAC Sigma 1 AC CLASS
VS3.12/C6	for	Sigma 1 AC CARDIO
VA1.1/C6	for	Sigma 1 AA STAR Sigma 1 AA CLASS
VA3.1/C6	for	Sigma 1 AA CARDIO

Sigma Doppler

Press the **SET** **MORE** **VERSION** touchkeys. After implementing the new software, the display will show :

V0.11/C6	for	Sigma Doppler
-----------------	-----	---------------

Proms exchange

Prom types and part numbers

Updated version : November 25th, 1988

Unit	Prom P.N.	Specification	Location	Type description
Sigma 1	866 903	630 082 - 12	Z 110	27128 K
	866 881	630 081 - 12	Z 109	27128 K
	866 911	630 083 - 12	Z 512	27128 K
	866 938	630 084 - 12	Z 513	27256
	866 946	630 085 - 12	Z 514	27256
Sigma 1 AC	866 903	630 082 - 12	Z 110	27128 K
	866 881	630 081 - 12	Z 109	27128 K
	866 911	630 083 - 12	Z 512	27128 K
	866 938	630 084 - 12	Z 513	27256
	866 946	630 085 - 12	Z 514	27256
Sigma 1 AA	875 236	630 118 - 2	Z 110	27128 K
	875 228	630 117 - 2	Z 109	27128 K
	875 244	630 119 - 2	Z 512	27256
	875 252	630 120 - 2	Z 513	27256
	875 279	630 121 - 2	Z 514	27256
Sigma Doppler	867 764	630 094 - 11	Z 110	27128 K
	867 772	630 095 - 11	Z 109	27128 K
	867 799	630 096 - 11	Z 512	27128 K
	867 802	630 097 - 11	Z 513	27128 K
	867 829	630 098 - 11	Z 514	27256 or 27C256
All units actual PAL's	866 954	630 086 - 1	Z 310	16 L 8
	858 714	630 027 - 1	Z 313	16 R 4
	866 962	630 087 - 1	Z 414	16 L 8
	865 117	630 056 - 3	Z 415	16 L 8

**Main improvements after implementing software version
V1.12/C6 or V3.12/C6 on Sigma 1**

1. Compatible with keyboard controller NEC-8279.
2. This software bears a lower servo loop gain at the control circuit for the wobbler motor.

**Main improvements after implementing software version
VS1.12/C6 or VS3.12/C6 on Sigma 1 AC**

1. Device configuration recognition implemented.
2. Compatible with keyboard controller NEC-8279.
3. This software bears a lower servo loop gain at the control circuit for the wobbler motor.
4. This software does not support any annular array transducers as defined for Sigma 1 AC.

**Main improvements after implementing software version
VA1.1/C6 or VA3.1/C6 on Sigma 1 AA**

1. Device configuration recognition implemented.
2. Some new parameters for some transducers.
3. Biopsy lines are implemented for all AA-transducers.
4. Compatible with keyboard controller NEC-8279.

**Main improvements after implementing software version
V0.11/C6 on Sigma Doppler**

1. Perivascular biometry is implemented "resistance Index and pulsality index".
2. Implementation of a wall filter of 800 Hz.

Product Group : Imaging
Ref. N° : Sigma Line - E.S.I. N° 5
Date : April 18th, 1983
Contact : Bernard MARTIN

Sigma Line

Engineering Service Information N° 5

Subject The "MONITOR OUT" signal at rear panel of Sigma 1 does not properly supply all the external video devices which are connected to Sigma .

Cause The "MONITOR OUT" of the Sigma may be slightly below the standard.

Solution Refer to "video output test procedure" page 2.

Modification If the voltage level is still too low despite the check of the video output as described on page 2 of 2, please contact Bernard Martin. Detailed information will be provided afterwards.

Caution Make sure that the interconnections between Sigma and the external video devices was performed according to the instructions given in the operating manual. Check that only 75 Ohm coaxial cables are used and that the line is loaded with one 75 Ohm device input only. All other inputs must be set to high impedance (Hi Z).

Additional information This modification should be implemented only if the described problem occurs and provided the conditions mentioned under "Caution" section are met.

Priority level is 3.

Units serial numbers 4010 and higher have been modified in the factory.

**Test of the
video output**

You need : to open the machine
a 3.5 MHz cardio probe
a 75 Ohm termination
an oscilloscope


Connect a 3.5 MHz cardio transducer

Set normal image
Post 2
TGC 1, TGC 2, ... , TGC 9 min
2D Gain min
TM gain min
Reject min

Observe the video signal of the Rec Out output terminated with a 75 Ohm resistor on the scope (200 mV/div, 10 us/div)


With the black porch as reference level, you should find the following voltage measurements within + 50 mV.

Overlay extra-white	710 mV
Image white for normal image	700 mV
Image white for inverse image	480 mV
Overlay white	480 mV
Image black	40 mV
Black porch	0 mV
Sync pulse	- 300 mV

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
Sigma 1 AA

Engineering Service Information N° 6

Subject Oscillations of the servoloop circuit on WOMOT p.c.b. when the annular array transducer is connected to Sigma 1 AA.
You can check this by connecting a scope on TP 12 (WOMOT p.c.b.) : a triangular oscillation of several volts at 20 to 25 kHz means that the p.c.b. must be modified.

Modification Insert a 220 pF 63 V ceramic capacitor at the place of C18 (see page 2 of 2).

Additional information Priority code: 2 (AA units should be modified at earliest convenience).
Units serial numbers 4128, 4134, 4135, 4141 and higher have been modified in the factory.
In order to avoid any problems when exchanging WOMOT p.c.b., we recommend that you immediately update your spares and your service case.
The modified WOMOT p.c.b. may be inserted into a Sigma 1 SC or a Sigma 1 without any problems.

Insert and solder C18 = 220pF (P.N. : 517 097)

Product Group : Imaging
Ref. N° : Sigma Line - E.S.I. N° 7
Date : April 18th, 1988
Contact : Bernard MARTIN

Sigma Line

Engineering Service Information N° 7

Subject When scanning with the 5 MHz Annular Array transducer in 10 cm depth, a black curved line appears at approximately 4 cm depth.

Cause The switching between two adjacent focussing areas is "visible" on the CRT in this particular configuration.

Solution Insert two zener diodes BZX55C3VO as described on page 2 of 2.

Caution This E.S.I. N° 7 applies for **Annular Array units only**, and not for Sigma 1 or Sigma 1 SC units.


Ordering information Diode BZX55C3VO
Part number : 732 753
Quantity : 2

Additional information

- . Priority code : 3 (the modification should be implemented on customer request only).
- . Sigma 1 AA serial numbers 4067, 4068, 4070, 4073, 4075, 4079, 4118 and higher have been modified in the factory.
- . If you do not feel confident to solder the two diodes on an SMD circuit, we will exchange the ATFOC f.o.c.

Modification

This is the new page 2 of 2 of Sigma Line
E.S.I. N° 7. Please take note that this page
cancels and replaces former page 2 of 2.

Product Group : Imaging
Instrument type/model : Sigma 1 Line (AA / AC)
Ref.: N° : Sigma 1 Line - E.S.I. N° 8
Date : November 6th, 1989
Contact : Eric CASTAING

SIGMA 1 LINE (AA / AC)

Engineering Service Information N° 8


Subject	New softwares release : SIGMA 1 AC LS (Class) - LSC (Class +) ____ VS 1.14/C6 SIGMA 1 AC CARDIO ____ VS 3.14/C6 SIGMA 1 AA LS - LSC ____ VA 1.3/C6 SIGMA 1 AA CARDIO ____ VA 3.3/C6
Improvements	<u>Transrectal Multiplane transducer compatibility :</u> - Stable at 0° / TM mode. - Right hand biopsy line - 105° / 211 beams image format. <u>Biopsy lines optimization :</u> - New biopsy line for AA 3.5 MHz and AA 5.0 MHz transducers. - Improvement of biopsy line for vaginal transducers. <u>Various improvements :</u> - Correction of "Distance/Ratio measurement in TM" bug, which used to lock the printer. - Correction of speed bug in TM.
Ordering information	One set of Proms will be shipped f.o.c. to each Kontron Company upon written request sent for the attention of Eric Castaing. Please, mention each needed version.
Additional information	Priority code is 3 : this modification should be implemented if a unit is at hand for repair or during a maintenance visit. You can check the E.S.I. level N° 8 on the label under the keyboard of your unit.

Kontron Instruments S.A.
2, Avenue de Manet - Montigny le Bretonneux
Boite Postale 81
78185 Saint Quentin en Yvelines Cedex
France


Tel : 33 (1) 30 43 81 52
Tlx : 698801 or 698802
Fax : 33 (1) 30 43 81 87

Softwares version check


Press the **SET** and **V** touchkeys and check that the display shows after installation :


For Sigma 1 AC LS - LSC

For Sigma 1 AC CARDIO

For Sigma 1 AA LS - LSC

For Sigma 1 AA CARDIO

Modification

1. Untighten the 6 locking screws on the base plate of Sigma 1 and disconnect the earth link between the frame and the top cover. Only then, remove the top cover.
2. Remove the metal cover of the left hand electronic box.
3. Remove all the ribbon and coaxial cables on the first three left hand pcb's.
4. Take SCON 2/MEMEX pcb's assy out of the box.
5. Exchange the proms on SCON 2 pcb (see locations on drawing page 4).
6. Perform re-assembly (from 5 to 1).
7. Check the unit performances.

Caution

This E.S.I. does not apply to SIGMA 1 SC.

Prom types and part numbers

Note : access time must be at least 200 ns.

VS 1.14/C6 software (set of 5 proms) for SIGMA 1 AC LS - LSC complete description				
Prom P.N.	Location	Specification	Checksum	Prom type
866 881	Z 109	630 081 -14	AA8D	27128
866 903	Z 110	630 082 -14	FBA0	27128
866 911	Z 512	630 083 -14	9C67	27128
866 938	Z 513	630 084 -14	8479	27256
866 946	Z 514	630 085 -14	8F3F	27256

VS 3.14/C6 software (set of 5 proms) for SIGMA 1 AC Cardio complete description				
Prom P.N.	Location	Specification	Checksum	Prom type
866 881	Z 109	630 081 -14	AA8D	27128
866 903	Z 110	630 082C -14	FBA1	27128
866 911	Z 512	630 083 -14	9C67	27128
866 938	Z 513	630 084 -14	8479	27256
866 946	Z 514	630 085 -14	8F3F	27256

VA 1.3/C6 software (set of 5 proms) for SIGMA 1 AA LS - LSC complete description				
Prom P.N.	Location	Specification	Checksum	Prom type
875 228	Z 109	630 117 - 3	AA8D	27128
875 236	Z 110	630 118 - 3	D51C	27128
875 244	Z 512	630 119 - 3	ACB4	27256
875 252	Z 513	630 120 - 3	48F2	27256
875 279	Z 514	630 121 - 3	2F1F	27256

VA 3.14/C6 software (set of 5 proms) for SIGMA 1 AA Cardio complete description				
Prom P.N.	Location	Specification	Checksum	Prom type
875 228	Z 109	630 117 - 3	AA8D	27128
875 236	Z 110	630 118C - 3	D51D	27128
875 244	Z 512	630 119 - 3	ACB4	27256
875 252	Z 513	630 120 - 3	48F2	27256
875 279	Z 514	630 121 - 3	2F1F	27256

SCON 2 pcb

SIGMA 1 AC

SP 630 081 - 14

27128

Z 109

SP 630 082 - 14

27128

Z 110

SP 630 083 - 14

27128

Z 512

SP 630 084 - 14

27256

Z 513

SP 630 085 - 14

27256

Z 514

SCON 2 pcb

SIGMA 1 AA

SP 630 117 - 3

27128

Z 109

SP 630 118 - 3

27128

Z 110

SP 630 119 - 3

27256

Z 512

SP 630 120 - 3


27256

Z 513

SP 630 121 - 3

27256

Z 514

Product Group : Imaging
Instrument type/model : Sigma 1 Line (AA / AC) US version
Ref.: N° : Sigma 1 Line - E.S.I. N° 10
Date : February 2nd, 1990
Contact : Eric CASTAING

SIGMA 1 LINE (AA / AC) US VERSION

Engineering Service Information N° 10

Subject

New softwares release :
SIGMA 1 AC LS (Class) - LSC (Class +) ____ VS 6.14/C6
SIGMA 1 AC CARDIO ____ VS 7.14/C6
SIGMA 1 AA LS - LSC ____ VA 6.3/C6
SIGMA 1 AA CARDIO ____ VA 7.3/C6

Improvements

Biometry modifications : 3 new tables
. HDC Head circumference
. ABC Abdominal circumference
. BND Binocular distance
Four former tables have been removed (THD ABD HIPA HIPB).

2 D to TM mode lock up problem is corrected.

Perfect horizontal tracing with the trackball.

New intrarectal multiplane probe frame rate : 8 4/7 NTSC

Ordering information

One set of Proms will be shipped f.o.c. to each Kontron Company upon written request sent for the attention of Eric Castaing. Please, mention each needed version.

Additional information

Priority code is **3** : this modification should be implemented if a unit is at hand for repair, during a maintenance visit or if a multiplane intrarectal probe is to be installed at customer's.

If the unit has been modified in the factory, the E.S.I. level N° 10 is mentioned on the "Révision/Modification" label under the keyboard of the unit.


If the modification is performed on the field, the service technician should tick off level N° 10 on this "Revision/Modification" label.

Kontron Instruments S.A.
2, Avenue de Manet - Montigny le Bretonneux
Boite Postale 81
78185 Saint Quentin en Yvelines Cedex
France


Tel. : 33 (1) 30 43 81 52
Tlx. : 698801 or 698802
Fax : 33 (1) 30 43 81 87

Softwares version check


Press the **SET** and **V** touchkeys and check that the display shows after installation :


For Sigma 1 AC LS - LSC

For Sigma 1 AC CARDIO

For Sigma 1 AA LS - LSC

For Sigma 1 AA CARDIO

Modification

Observe electrostatic discharge safety rules during these operations.

1. Untighten the 6 locking screws on the base plate of Sigma 1 and disconnect the earth link between the frame and the top cover. Only then, remove the top cover.
2. Remove the metal cover of the left hand electronic box.
3. Remove all the ribbon and coaxial cables on the first three left hand pcb's.
4. Take SCON 2/MEMEX pcb's assy out of the box.
5. Exchange the proms on SCON 2 pcb (see locations on drawing page 4).
6. Perform re-assembly (from 5 to 1).
7. Check the unit performances.
8. Tick off N° 10 on E.I.S level label (Révision/Modification) under the keyboard.
9. Check that the customer's Operating Manual comes within this modification.

Caution

This E.S.I. does not apply to SIGMA 1 SC.

Prom types and part numbers

Note : access time must be at least 200 ns.

VS 6.14/C6 software (set of 5 proms) for SIGMA 1 AC LS - LSC complete description

Prom P.N.	Location	Specification	Checksum	Prom type
888 881	Z 109	630 081 -14	1CAA8D	27128
415 243	Z 110	630 310 - 1	1CFD20	27128
415 251	Z 512	630 311 - 1	1BA7E7	27128
415 278	Z 513	630 312 - 1	41881A	27258
415 288	Z 514	630 313 - 1	34FE0C	27258

VS 7.14/C6 software (set of 5 proms) for SIGMA 1 AC Cardio complete description

Prom P.N.	Location	Specification	Checksum	Prom type
888 881	Z 109	630 081 -14	1CAA8D	27128
415 332	Z 110	630 318 - 1	1CFD20	27128
415 251	Z 512	630 311 - 1	1BA7E7	27128
415 278	Z 513	630 312 - 1	41881A	27258
415 288	Z 514	630 313 - 1	34FE0C	27258

VA 6.3/C6 software (set of 5 proms) for SIGMA 1 AA LS - LSC complete description






Prom P.N.	Location	Specification	Checksum	Prom type
875 228	Z 109	630 117 - 3	1CD51C	27128
415 318	Z 110	630 318 - 1	1CDCA2	27128
415 294	Z 512	630 314 - 1	21B8F0	27258
415 308	Z 513	630 315 - 1	396AE9	27258
415 324	Z 514	630 317 - 1	35A0C8	27258

VA 7.14/C6 software (set of 5 proms) for SIGMA 1 AA Cardio complete description

Prom P.N.	Location	Specification	Checksum	Prom type
875 228	Z 109	630 117 - 3	1CD51C	27128
415 340	Z 110	630 319 - 1	1CDCA2	27128
415 294	Z 512	630 314 - 1	21B8F0	27258
415 308	Z 513	630 315 - 1	396AE9	27258
415 324	Z 514	630 317 - 1	35A0C8	27258






SCON 2 pcb


SIGMA 1 AC

SCON 2 pcb

SIGMA 1 AA

SIGMA 1

Engineering Service Information N° 11

Subject NEC cathod ray tube of the L&B monitor (P.N. 861 561) is not manufactured any more.

Solution Use new L&B monitor with CLINTON cathod ray tube.

Ordering information	P.N. 898 430	L&B Clinton monitor kit (Sigma 1 E.S.I. N° 11) including :
		1 x L&B Clinton monitor 410 403
		4 x Washes DD3 407 712
		2 x CRT holder 408 093

Additional information Priority code is 4 (information). As the TVM pcb is matched to the CRT, it is necessary to replace the whole assembly.


Units serial numbers 5270 and higher have been modified in the factory.

Caution Observe electrostatic discharge protection rules.

Kontron Instruments S.A.
2, Avenue du Manet - Montigny le Bretonneux
Boite Postale 81
78185 Saint Quentin en Yvelines Cedex
France

Tel : 33 (1) 30 43 81 52
Tlx : 698801 or 698802
Fax : 33 (1) 30 44 23 57

Engineering Service Information

SIGMA 1 AC/AA (all versions)

Engineering Service Information N° 12

Subject	1.	Noise on ECG signal.
	2.	Black lines radiating on 2D picture.
Solution	Hardware modification (see modification page 2) :	
	1.	Resistor 1 KOhm
	2.	Resistor 330 Ohm
Ordering information	1.	Resistor 1 KOhm 5 % 1/4 W P.N. 510 408
	2.	Resistor 330 Ohm 5 % 1/4 W P.N. 510 343
Additional information	Priority code is 3. This modification should be implemented if a unit is at hand for repair, during a maintenance visit, or if the problem is raised by a customer.	
	Units serial numbers 5369 and higher have been modified in the factory.	
Caution	Observe electrostatic discharge protection rules.	

Kontron Instruments S.A.
2, Avenue du Manet - Montigny le Bretonneux
Boite Postale 81
78185 Saint Quentin en Yvelines Cedex
France

Tel.: 33 (1) 30 43 81 52
Tlx. : 696801 or 696802
Fax: 33 (1) 30 44 23 57

Engineering Service Information

Modification

1. E.C.G. module

Remove ECG module from SIGMA 1.

Connect a 1 K resistor between IC3 pin 14 (+5V) and pins 8 and 9 (gate inputs) as a pull up resistor.

Reassemble in reverse order.

Perform safety checks.

Amend your service manual.

2. AAFRO p.c.a.

Remove AAFRO p.c.a.

Solder a 330 Ohm resistor between Z 409 / pin 23 and + 5V rails as indicated on the schematics and drawings.

Reassemble in reversed order.

Perform safety and performance checks.

Amend your service manual.

Product Group: **Imaging**
Ref. No: **90511 K.I.S. EC/mm**
Date: **July 6th, 1990**
Contact: **Eric CASTAING**

SIGMA 1 LINE - E.C.G. MODULE

Dear colleagues,

The former ECG module P.N. 855 413 is replaced by ECG module P.N. 889 687 which is the one we use for SIGMA 44 family.

The difference lies in the safety PSU of the ECG p.c.b. which is replaced by a piggy back PSU module and a pull-up resistor as explained in E.S.I. SIGMA 1 N° 12.

This new module will prevent the equipment from current spikes on + 15 V rails and from noise on ECG signal.

SIGMA 1 units serial number 5379 and higher have been modified in the factory.

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


Product Group: **Imaging**
 Ref. No: **90388 K.I.S. EC/mm**
 Date: **June 13th, 1990 (second issue)**
 Contact: **Eric CASTAING**

SIGMA 1 - PROVISIONNAL SOFTWARE

It happens that for some special cases, we may deliver accessories with " provisional " softwares. This means that Eproms could be labeled under a " software module " number instead of a specification number.

For your information, here is the location on SCON 2 p.c.b.

SIGMA 1 LINE

Engineering Service Information N° 13

Subject New sets of German software. The first available versions are :
VS 8.14 VS 9.14 VA 8.3 VA 9.3

New features Introduction of Dr. Hansmann OB/GYN biometry tables in SIGMA with dedicated software sets. Former OB/GYN tables are removed.

Ordering information Refer to E.S.I. N° 14..

Additional information From October 1990, all SIGMA 1 software releases will include this German version as explained in next SIGMA 1 LINE E.S.I. N° 14.

Tick off N° 13 on status label of either versions 9.xx or 8.xx.

Caution This E.S.I. does not apply to SIGMA 1 SC.

Product Group: Imaging
Ref. No: 90997 K.I.S. EC/mm
Date: December 3rd, 1990
Contact: Eric CASTAING

SIGMA 1 : E.S.I. N° 14

Since the release of this E.S.I., the eprom 630 081 has been modified on SIGMA 1 AC software. It is now version 630 081-15.


The software revision number remains unchanged.

Please, make sure you have this eprom exchanged on the following units : 5636 5637 5643 5644 5647 5648 5649 5641 5652 5653 5654 5655 5656 5657 5658 5659 5660 5661 5662 5663 5664 5665 5671 5688.

This change is notified in SIGMA 1 E.S.I. N° 14 and in SIGMA DOPPLER E.S.I. N° 15 dated November 22nd, 1990 (second edition).

With best regards,

Eric CASTAING

p/d

Page 1/1

Kontron Instruments S.A.,
2 Avenue de Manet,
Montigny-Le Bretonneux,
(Yvelines),
France
Telephone 03-0438150
Telex 303115

Kontron Instruments S.p.A.,
Via G. Fantoli 16/15
20138 Milano
Italy
Telephone 02-50721
Telex 312288

Kontron Instruments A.G.,
Bernstr. Sud 169,
8010 Zurich,
Switzerland.
Telephone 01-4354111
Telex 522191 A

Kontron Instruments Ltd.,
Blackmoor Lane
Craxley Centre,
Watford, Herts. WD1 8XQ,
United Kingdom
Telephone 0923 245991/56252
Telex 332712

SIGMA 1 LINE

Engineering Service Information N° 14

Subject New software release.

New features

7.5 S transducer optimization.

2D to TM mode lock up problem is corrected (*).

Perfect horizontal tracing with trackball (*).

New multiplane intra-rectal probe frame rate (*) from 10 fps down to 8 1/3 PAL (84/7 NTSC).

Random inverse video display in imaging mode (WGA gate array bug) corrected.

(*) already corrected in US version. Refer to E.S.I. N° 10.

Ordering information Each country will be provided with one set of required software version(s) upon written request sent for the attention of Lilliane BESNARD.

Kit soft AC/AC Cardio standard version	888 230
Kit soft AA/AA Cardio standard version	888 540

Kit soft AC/AC Cardio USA	version	888 370
Kit soft AA/AA Cardio USA	version	888 680

Kit soft AC/AC Cardio German	version	898 740
Kit soft AA/AA Cardio German	version	898 910

Additional information Priority code is 3 . Units S.N. 5636 and higher have been modified in the factory. This E.S.I. does not apply to SIGMA 1 SC.

Kontron Instruments S.A.
2, Avenue du Mahet - Montigny le Bretonneux
Boite Postale 81
78185 Saint Quentin en Yvelines Cedex
France

Tel. : 33 (1) 30 43 81 52
Tlx. : 698801 or 698802
Fax : 33 (1) 30 43 81 87

Engineering Service Information

Modification

Observe electrostatic discharge safety rules during these operations.

1. Untighten the 6 locking screws on the base plate of SIGMA 1 and disconnect the earth link between the frame and the top cover. Only then, remove the top cover.
2. Remove the metal cover of the left hand electronic box.
3. Remove the ribbon and coaxial cables on the first three left hand pcb's.
4. Take SCON 2/ MEMEX pcb's assy out of the box.
5. Exchange the prompts on SCON 2 pcb (see locations on drawing pages 3 to 5).
6. Perform re-assembly (from 5 to 1).
7. Check the unit performances.
8. Tick off N° 14 on E.S.I. level label (Revision/Modification) under the keyboard.
9. Tick off N° 13 on E.S.I. level label (Revision/Modification) under the keyboard if you use German software (Hansmann).
10. Check that the customer's Operating Manual comes within this modification.

Software version check

Press the  and  touchkeys and check that the screen shows one of the following displays after installation :

VS 1.15/C6
SIGMA 1 AC LS - LSC (Standard)

VS 3.15/C6
SIGMA 1 AC CARDIO (Standard)

VA 1.4/C6
SIGMA 1 AA LS - LSC (Standard)

VA 3.4/C6
SIGMA 1 AA CARDIO (Standard)

VS 6.15/C6
SIGMA 1 AC LS - LSC (US)

VA 7.15/C6
SIGMA 1 AA CARDIO (US)

VA 6.4/C6
SIGMA 1 AA LS - LSC (US)

VA 7.4/C6
SIGMA 1 AA CARDIO (US)

VS 8.15/C6
SIGMA 1 AC LS - LSC (German)

VS 9.15/C6
SIGMA 1 AA CARDIO (German)

VA 8.4/C6
SIGMA 1 AC LS - LSC (German)

VA 9.15/C6
SIGMA 1 AA CARDIO (German)

Prom types and part numbers of STANDARD version

Note : access time must be a least 200 ns.

VS 1.15/C6	software (set of 5 proms) for SIGMA 1 AC				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	888 881	Z 109	630 081-14	00ICAA8D	27128
	888 903	Z 110	630 082-15	00ICFF68	27128
	888 911	Z 512	630 083-15	001BA33D	27128
	888 938	Z 513	630 084-15	0042084E	27256
	888 945	Z 514	630 085-15	0034C0E7	25256

VS 3.15/C6	software (set of 5 proms) for SIGMA 1 AC CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	888 881	Z 109	630 081-15	00ICAA8D	27128
	419 052	Z 110	630 353-2	00ICFF67	27128
	888 911	Z 512	630 083-15	001BA33D	27128
	888 938	Z 513	630 084-15	00420E4E	27256
	888 945	Z 514	630 085-15	0034C0E7	25256

VA 1.4/C6	software (set of 5 proms) for SIGMA 1 AA				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	630 117-4	00IC8805	27128
	875 236	Z 110	630 118-4	00ICDCC4	27128
	875 244	Z 512	630 119-4	0021AA5F	27256
	875 252	Z 513	630 120-4	0039DFF8	27256
	875 279	Z 514	630 121-4	00352CD0	25256

VA 3.4/C6	software (set of 5 proms) for SIGMA 1 AA CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	630 117-4	00IC8805	27128
	419 745	Z 110	630 380	001CDCC5	27128
	875 244	Z 512	630 119-4	0021AA5F	27128
	875 252	Z 513	630 120-4	0039DFF8	27256
	875 279	Z 514	630 121-4	00352CD0	25256

Prom types and part numbers of US version

Note : access time must be a least 200 ns.

VS 6.15/C6	software (set of 5 proms) for SIGMA 1 AC				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	888 881	Z 109	630 081-14	00ICAA8D	27128
	415 243	Z 110	630 310-2	001CFF6B	27128
	415 251	Z 512	630 311-2	001BA33D	27128
	415 278	Z 513	630 312-2	00420649	27256
	415 288	Z 514	630 313-2	0035273C	25256

VS 7.15/C6	software (set of 5 proms) for SIGMA 1 AC CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	888 881	Z 109	630 081-14	00ICAA8D	27128
	415 332	Z 110	630 318-2	001CFF6B	27128
	415 251	Z 512	630 311-2	001BA33D	27128
	415 278	Z 513	630 312-2	00420649	27256
	415 288	Z 514	630 313-2	0035273C	25256

VA 6.4/C6	software (set of 5 proms) for SIGMA 1 AA				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	630 117-4	00IC8805	27128
	415 318	Z 110	630 316-2	001CDCC9	27128
	415 294	Z 512	630 314-2	0021AA5F	27256
	415 308	Z 513	630 315-2	0039E0F8	27256
	415 324	Z 514	630 317-2	003592E8	25256

VA 7.4/C6	software (set of 5 proms) for SIGMA 1 AC CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	630 117-4	00IC8805	27128
	415 340	Z 110	630 319-2	001CDCC9	27128
	415 294	Z 512	630 314-2	0021AA5F	27256
	415 308	Z 513	630 315-2	0039E0F8	27256
	415 324	Z 514	630 317-2	003592E8	25256

Prom types and part numbers of GERMAN version

Note : access time must be a least 200 ns.

VS 8.15/C6	software (set of 5 proms) for SIGMA 1 AC				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	886 881	Z 109	830 081-14	001CAA8D	27128
	419 109	Z 110	830 354-2	001CFF8D	27128
	419 117	Z 512	830 355-2	001BA33D	27128
	419 141	Z 513	830 356-2	004218EA	27258
	419 168	Z 514	830 357-2	0034C733	25258


VS 9.15/C6	software (set of 5 proms) for SIGMA 1 AC CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	886 881	Z 109	830 081-14	001CAA8D	27128
	419 178	Z 110	830 358-2	001CFF8D	27128
	419 117	Z 512	830 355-2	001BA33D	27128
	419 141	Z 513	830 356-2	004218EA	27258
	419 168	Z 514	830 357-2	0034C733	25258

VA 8.4/C6	software (set of 5 proms) for SIGMA 1 AA				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	830 117-4	001C8805	27128
	419 184	Z 110	830 359-2	001C0CC8	27128
	419 192	Z 512	830 360-2	0021AA5F	27258
	419 208	Z 513	830 361-2	0039DFFE	27258
	419 214	Z 514	830 362-2	0035330F	25258

VA 9.4/C6	software (set of 5 proms) for SIGMA 1 AC CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	830 117-14	001C8805	27128
	419 222	Z 110	830 363-2	001C0CC8	27128
	419 192	Z 512	830 360-2	0021AA5F	27258
	419 208	Z 513	830 361-2	0039DFFE	27258
	419 214	Z 514	830 362-2	0035330F	25258


SCON 2 pcb


SIGMA 1 AC

SCON 2 pcb

SIGMA 1 AA

Product Group : Imaging
Instrument type/model : SIGMA 1 LINE
Ref.: N° : SIGMA DOPPLER E.S.I. N° 15
Date : November 22nd, 1990
Contact : Eric CASTAING

SIGMA 1 LINE - SIGMA DOPPLER

Engineering Service Information N° 15

CANCELS AND REPLACES EDITION ISSUED ON OCTOBER 18th, 1990

Subject

New software release. Also see E.S.I. SIGMA DOPPLER N° 15.

New features

Optimization for 7.5 S transducer.

Modification

Install the new software on DOPPLER SCON 2 pca.

Check the display shows **V0.12/C6** after implementation of the new software
by pressing the **SET** **MORE** **VERSION** touchkeys.

Ordering information

Each country will be provided with one set of required software version(s) upon
written request sent for the attention of Liliane BESNARD.

Kit soft SCON2 Sigma Doppler : P.N. 413 380

Additional information

Priority code is 3 . Units S.N. 880 and higher have been modified in the factory.


Compatibility level : C6.

Prom types and part numbers of V0.12/C6 version

Note : access time must be a least 200 ns.

V0.12/C6	software (set of 5 proms) for SIGMA DOPPLER			complete description	
	Prom P.N.	Location	Specification	Checksum	Prom type
	867 764	Z 110	630 094-12	002D9EC5	27128
	867 772	Z 109	630 093-11	003ADE62	27128
	867 799	Z 512	630 095-12	001B551A	27128
	867 802	Z 513	630 096-12	000DFCC5	27128
	867 829	Z 514	630 097-12	003AC31D	25258

Proms location on Doppler SCON 2 pca

CANCELS AND REPLACES EDITION ISSUED ON OCTOBER 18th, 1990

SIGMA 1 LINE

Engineering Service Information N° 16

Subject

New software release.

New features

New multiplane endorectal probe (Capistrano) P.N. 421 480 optimization.

Ordering information

Each country will be provided with one set of required software version(s) upon written request sent for the attention of Liliane BESNARD.

Kit soft AC/AC Cardio standard	version	888 230
Kit soft AA/AA Cardio standard	version	888 540

Kit soft AC/AC Cardio USA	version	888 370
Kit soft AA/AA Cardio USA	version	888 680

Kit soft AC/AC Cardio German	version	898 740
Kit soft AA/AA Cardio German	version	898 910

Additional information

Priority code is 3. Units S.N. 6139, 6140, 6143, 6144, 6147 and higher have been modified in the factory. Compatibility level = C6. This E.S.I. does not apply to SIGMA 1 SC.

Observe electrostatic protection rules when handling subassemblies !

Modification

Observe electrostatic discharge safety rules during these operations.

1. Untighten the 6 locking screws on the base plate of SIGMA 1 and disconnect the earth link between the frame and the top cover. Only then, remove the top cover.
2. Remove the metal cover of the left hand electronic box.
3. Remove the ribbon and coaxial cables on the first three left hand pcb's.
4. Take SCON 2/ MEMEX pcb's assy out of the box.
5. Exchange the proms on SCON 2 pcb (see locations on drawing pages 3 to 5).
6. Perform re-assembly (from 5 to 1).
7. Check the unit performances.
8. Tick off N° 16 on E.S.I. level label (Revision/Modification) under the keyboard.
9. Tick off N° 13 on E.S.I. level label (Revision/Modification) under the keyboard if you use German software (Hansmann).
10. Check that the customer's Operating Manual comes within this modification.

Software version check

Press the **SET** and **V** touchkeys and check that the screen shows one of the following displays after installation :

VS 1.16/C6

SIGMA 1 AC LS - LSC (Standard)

VS 3.16/C6

SIGMA 1 AC CARDIO (Standard)

VA 1.5/C6

SIGMA 1 AA LS - LSC (Standard)

VA 3.5/C6

SIGMA 1 AA CARDIO (Standard)

VS 6.16/C6

SIGMA 1 AC LS - LSC (US)

VS 7.16/C6

SIGMA 1 AC CARDIO (US)

VA 6.5/C6

SIGMA 1 AA LS - LSC (US)

VA 7.5/C6

SIGMA 1 AA CARDIO (US)

VS 8.16/C6

SIGMA 1 AC LS - LSC (German)

VS 9.16/C6

SIGMA 1 AC CARDIO (German)

VA 8.5/C6

SIGMA 1 AA LS - LSC (German)

VA 9.5/C6

SIGMA 1 AA CARDIO (German)

Prom types and part numbers of STANDARD version

Note : access time must be a least 200 ns.

VS 1.16/C6	software (set of 5 proms) for SIGMA 1 AC				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	868 881	Z 109	830 081-15	001C24E5	27128
	868 903	Z 110	830 082-17	001C850A	27128
	868 911	Z 512	830 083-17	001B9F8D	27128
	868 938	Z 513	830 084-17	0042011D	27258
	868 945	Z 514	830 085-17	0034AC07	25258
VS 3.16/C6	software (set of 5 proms) for SIGMA 1 AC CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	868 881	Z 109	830 081-15	001C24E5	27128
	419 052	Z 110	830 353-4	001C850B	27128
	868 911	Z 512	830 083-17	001B9F8D	27128
	868 938	Z 513	830 084-17	0042011D	27258
	868 945	Z 514	830 085-17	0034AC07	25258
VA 1.5/C6	software (set of 5 proms) for SIGMA 1 AA				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	830 117-4	001C8805	27128
	875 238	Z 110	830 118-8	001C8EL82	27128
	875 244	Z 512	830 119-8	0021C077	27258
	875 252	Z 513	830 120-8	00387C48	27258
	875 279	Z 514	830 121-8	00351731	25258
VA 3.5/C6	software (set of 5 proms) for SIGMA 1 AA CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	830 117-4	001C8805	27128
	419 745	Z 110	830 380-3	001C8EL82	27128
	875 244	Z 512	830 119-8	0021C077	27258
	875 252	Z 513	830 120-8	00387C48	27258
	875 279	Z 514	830 121-8	00351731	25258

Prom types and part numbers of US version

Note : access time must be a least 200 ns.

VS 6.16/C6	software (set of 5 proms) for SIGMA 1 AC				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	868 881	Z 109	830 081-15	001C24E5	27128
	415 243	Z 110	830 310-4	001C850F	27128
	415 251	Z 512	830 311-4	001B9F8D	27128
	415 278	Z 513	830 312-4	0042011D	27258
	415 288	Z 514	830 313-4	00354A07	25258
VS 7.16/C6	software (set of 5 proms) for SIGMA 1 AC CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	868 881	Z 109	830 081-15	001C24E5	27128
	415 332	Z 110	830 318-4	001C850F	27128
	415 251	Z 512	830 311-4	001B9F8D	27128
	415 278	Z 513	830 312-4	0042004D	27258
	415 288	Z 514	830 313-4	00351245	25258
VA 6.5/C6	software (set of 5 proms) for SIGMA 1 AA				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	830 117-4	001C8805	27128
	415 318	Z 110	830 318-4	001C8E87	27128
	415 294	Z 512	830 314-4	0021C07A	27258
	415 308	Z 513	830 315-4	00387C48	27258
	415 324	Z 514	830 317-4	00357C3F	25258
VA 7.5/C6	software (set of 5 proms) for SIGMA 1 AA CARDIO				complete description
	Prom P.N.	Location	Specification	Checksum	Prom type
	875 228	Z 109	830 117-4	001C8805	27128
	415 340	Z 110	830 319-4	001C8E87	27128
	415 294	Z 512	830 314-4	0021C07A	27258
	415 308	Z 513	830 315-4	00387C48	27258
	415 324	Z 514	830 317-4	00357C3F	25258

Prom types and part numbers of GERMAN version

Note : access time must be a least 200 ns.

VS 8.16/C8 software (set of 5 proms) for SIGMA 1 AC complete description

Prom P.N.	Location	Specification	Checksum	Prom type
868 881	Z 109	830 081-15	001C24E5	27128
419 109	Z 110	830 354-4	001CB511	27128
419 117	Z 512	830 355-4	001B58FD	27128
419 141	Z 513	830 356-4	0042011D	27258
419 188	Z 514	830 357-4	0034B248	25258

VS 9.16/C8 software (set of 5 proms) for SIGMA 1 AC CARDIO complete description

Prom P.N.	Location	Specification	Checksum	Prom type
868 881	Z 109	830 081-15	001C24E5	27128
419 176	Z 110	830 358-4	001CB511	27128
419 117	Z 512	830 355-3	001B58FD	27128
419 141	Z 513	830 356-3	0042011D	27258
419 188	Z 514	830 357-3	0034B248	25258

VA 8.5/C8 software (set of 5 proms) for SIGMA 1 AA complete description


Prom P.N.	Location	Specification	Checksum	Prom type
875 228	Z 109	830 117-4	001C8805	27128
419 184	Z 110	830 358-4	001C8E59	27128
419 192	Z 512	830 380-4	0021C07A	27258
419 206	Z 513	830 381-4	00397C48	27258
419 214	Z 514	830 382-3	00351D72	25258

VA 9.5/C8 software (set of 5 proms) for SIGMA 1 AA CARDIO complete description

Prom P.N.	Location	Specification	Checksum	Prom type
875 228	Z 109	830 117-4	001C8805	27128
419 222	Z 110	830 383-4	001C8E59	27128
419 192	Z 512	830 380-3	0021C07A	27258
419 206	Z 513	830 381-3	00397C48	27258
419 214	Z 514	830 382-3	00351D72	25258


SCON 2 pcb


SIGMA 1 AC

SCON 2 pcb

SIGMA 1 AA

Product Group: Imaging
Ref. No: 90314 K.I.S. EC/mm
Date: September 14th, 1990
Contact: Eric CASTAING

MODIFICATION ON PERIP P.C.B.

The first version of PERIP p.c.b. does not include the ECG module detection nor the CARDIS signal inhibition which may cause some spurious problem with Doppler.

A. If you wish to upgrade this p.c.b. on site, here is the procedure :
(schematics ref. 03 121 006 and 03 121 B 306)

1. Disconnect J 811 / Pin 15 from ground on both sides of the p.c.b.
2. Link J 811 / Pin 15 to C21 and R 30.
3. Cut the tracks from C 21 - R 30 to J 810 / Pin 34 (near the connector).
4. Connect 1 resistor 10 KOhm 1/4 W 5 % between Z 105 / Pin 11 to + 5V.
Connect 1 resistor 10 KOhm 1/4 W 5 % between Z 105 / Pin 13 to + 5V.

All old p.c.b.'s returned to the factory are upgraded.


B. If you need to connect a SIGMA DOPPLER to a SIGMA 1 which is not a Cardio version, then you have to simulate the ECG module to enable the DOPPLER operation. The modification is as follows :

Link J811 / PIN 15 to 0 Volt.

With best regards,

Eric CASTAING

p/p

Page 1/1

Kontron Instruments S.A.,
2 Avenue de Manet
Montigny-Le Bretonneux
(Yvelines)
France
Telephone 03-0436152
Telex 890500

Kontron Instrumente S.p.A.,
Via G. Fantoli 16/15
20138 Milano
Italy
Telephone 03-50721
Telex 312288

Kontron Instruments A.G.,
Bernerstr. Sud 169,
8010 Zurich,
Switzerland
Telephone 01-4354111
Telex 822191 A

Kontron Instruments Ltd.,
Blackmoor Lane
Croxley Centre
Watford, Herts. WD1 8XO
United Kingdom
Telephone 0920 245991 56362
Telex 930701

Product Group: Imaging
Ref. No: 90586 K.I.S. EC/mm
Date: September 14th, 1990
Contact: Eric CASTAING

ULTRASOUND TRANSDUCERS : IMPROVEMENTS

Dear Colleagues,

Since June, all mechanical scanheads produced in Montigny factory have been improved with a stainless steel ring around the membrane in the contact area with the body. This new ring increases the tightness and consequently reduces the occurrence of unstuck membrane and bubbles.


Transducers serial numbers listed herebelow and higher have been modified in the factory :

Wobbler	3.5	A	# 1858	Annular Array	3.5	A	# 10342
Wobbler	3.5	B	# 3768	Annular Array	3.5	B (D)	# 18203
Wobbler	3.5	C	# 33182	Annular Array	5.0	A	# 13263
Wobbler	3.5	D	# 7437	Annular Array	5.0	B (D)	# 19125
Wobbler	3.5	DS	# ALL	Annular Array	7.5	B (D)	# 20141
Wobbler	7.5	B	# 5632				

All transducer sent out for repair will be upgraded.

With best regards,

Eric CASTAING

P/q

Page 1/1

Kontron Instruments S.A.,
2 Avenue de Manet,
Montigny Bretonneux,
(Yvelines)
France
Telephone 01-69-68-152
Telex 70663

Kontron Instruments S.p.A.,
Via G. Fantini 16/15
20138 Milano
Italy
Telephone 02-50701
Telex 212235

Kontron Instruments A.G.,
Bernstr. Süd 169,
8010 Zurich,
Switzerland
Telephone 01-4354111
Telex 822191 A

Kontron Instruments Ltd.,
Blackmoor Lane
Croxley Centre
Watford, Herts. WD1 8XQ,
United Kingdom
Telephone 0923 245991/56352
Telex 332212

Product Group: Imaging
Ref. No: 90663 K.I.S. EC/mm
Date: September 14th, 1990
Contact: Eric CASTAING

SIGMA 1 KEYBOARDS

This is to inform you that the "flat" keyboard is no longer available. Only the "tactile" keyboard remains.

Tactile keyboard English	P.N. 874 523
Tactile keyboard French	P.N. 874 604
Keyboard ribbon cable 18B	P.N. 884 413
Connector shell	P.N. 884 421

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

Spark
Nut
Electrolytic
Collar
Coated/Covered
Input(s)
And
Tin
Label
Extractor

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Spacer
Interface
Interrupter
Inverter
Insulation
Insulated

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Wire carrier
 Peripheral
 Cell/Battery
 Flat
 Stud
 Plastic
 Holder
 Potentiometer
 Push button
 Plug/Socket
 Programmable
 Shelding
 PVC

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RAM
Reference
Regulator
Relay
Alternative replacement
Guiding mark
Network
Resistor
Shrink
Retriggerable
Rivet
Round
Washer
Red

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Shotky
Selector
Sequensor
Locked
To lock
SIL
Simple/Single/Mono
Solder
Soldered
Flexible
turn
Socket

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